# Invacare® Kompas®

- en Power Wheelchair User Manual**
- fr Fauteuil roulant électrique Manuel d'utilisation**
- da El-køretøj Brugsanvisning**
- it Carrozzina elettrica Manuale d'uso**
- es Silla de ruedas eléctrica Manual del usuario**
- nl Elektrische rolstoel Gebruikershandleiding**
- pt Cadeira de rodas elétrica Manual de Utilização**
- de Elektrorollstuhl Gebrauchsanweisung**

This manual MUST be given to the user of the product.  
BEFORE using this product, this manual MUST be read and saved for  
future reference.

**Yes, you can.®**

© 2023 Invacare Corporation

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by™ and®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

# User Manual

## Table of Content

<b>Introduction</b>	
Practice before operating	2
<b>Safety Instructions</b>	
Considerations	4
Driving Outside	5
Driving on various terrains	7
<b>Familiarize yourself with your power wheelchair</b>	
Feature Diagram	9
General Specification Data	10
Assembly	11
<b>Adjustment for comfort</b>	
Adjustment of Joystick	14
<b>Operation</b>	
Controls and Indicators	15
Driving, Steering and Braking	15
Pushing the Powerchair by Hand	17
Shutdown	18
Getting Ready for Permanent use	18
<b>Batteries</b>	19
<b>Battery Charging</b>	21
<b>Lithium Battery Instructions for Use</b>	22
<b>Folding the Wheelchair</b>	24
<b>Import and Safety Information</b>	25
P&G VR2 Controller	28
<b>Maintenance</b>	
Rear Wheel	33
Front Wheel	33
<b>Troubleshooting Guide</b>	34
<b>IEC Symbols</b>	35
<b>Disclosure Information(ISO)</b>	36
<b>Warranty</b>	37

## **Introduction**

Welcome aboard your new Kompas Power Wheelchair and thank you for choosing our product. This latest model has been designed with specific practical user needs in mind. It combines solid rugged construction and modern high-tech electronics, safety and performance.

The Kompas Power Wheelchair has a programmable electronic control system reflecting the latest state of the art design. It means, among other things, that your Power Wheelchair can be programmed and adjusted within a given range of its performance characteristics, to suit your individual needs and wishes. The controller is set at the factory to give the Power Wheelchair nominal operating performance characteristics.

If after an initial trial period, when you have familiarized yourself with your powerchair, you would prefer certain features to be more adapted to your personal requirements (e.g. acceleration, deceleration, maximum speed, cornering speed, response to joystick movement ,etc.) Please contact your local dealer. They will be able to advise you, and they have the necessary equipment to carry out any adjustments required.

We recommend that your power wheelchair is serviced by your local dealer every six months. It is in your best interests, not only to ensure your personal safety but also to ensure long life and reliability. Please ensure that your service schedule (at the rear of this booklet) is signed by your dealer after each service.

This manual contains some useful tips and information on safety, operation and maintenance. Please read it carefully to ensure that you get maximum enjoyment and benefit from your new independence and mobility. Users with visual impairments can have the documentation read out by a helper.

Whenever you require special advice and attention, do not hesitate to contact your local dealer. They have all the equipment and know-how to provide expert service. You are their customer as much as ours and they want to see you satisfied.

In case of a serious incident with the product, you should inform the manufacturer and the competent authority in your country.

# User Manual


**Failure to follow these instructions may result in damage to the power wheelchair or serious injury.**

## ■ Practice Before Operating


Find an open area such as a park and have an assistant to help you practice until you have confidence operating this vehicle.

Make sure that the power is off before getting in or out of the seat. Set the speed control button according to your driving ability.

**We recommend that you keep the speed control at the slowest position until you are familiar with the driving characteristics of this vehicle. Users should have no hearing and vision impairment, and responsive.**

## Getting familiar with this vehicle

**First, practice moving forward.  
Be sure to set the speed to the lowest setting.**

**After becoming familiar with moving forward, practice marking "S" turns.**

**Once you are familiar with "S" turns, practice moving in reverse. Note that for any speed control setting, the vehicle moves more slowly in reverse than forward.**

# User Manual

## ■ Safety Considerations

DO NOT do any of the following

**NO!**

Do not carry any passengers

**NO!**

Do not drive across a slope

**NO!**

Do not drink and drive  
Consult your physician to  
determine if your medications  
impair your ability to control this  
vehicle

**NO!**

Do not tow a trailer

**NO!**

Do not turn on or use hand-held personal  
communication devices such as citizens band(CB)  
radios and cellular phones

## ■ Driving Outdoors

When you are on the road, please pay attention to the following:

**NO!**

Do not drive in traffic.

**NO!**

Do not drive beside a river, port, or lake without a fence or railing.

**NO!**

If possible, do not drive during the rain.

**NO!**

If possible, do not drive during or on snow.

**NO!**


Do not drive off-road or on any uneven surfaced roads.

**NO!**

If possible, do not drive at night.

# User Manual

**NO!**

Make sure that there are no obstacles behind you when in reverse.

We recommend to set up the speed at the lowest setting for reversing.

**NO!**

Do not make sudden stops, weave erratically, or make sharp turns.

**NO!**

Keep your arms on or inside the armrests and feet on the footrest at all time.

**NO!**

Do not attempt to climb curbs greater than 1.2"(3cm).

## ■ Driving on Various Terrains

Driving on hills is more dangerous than on level surfaces. If you fail to heed these warnings, a fall, tip-over or loss of control may occur and cause severe injury to the vehicle user or others.

**NO!**

**Do not attempt to climb a hill greater than 6°**

**NO!**

**Do not reverse while driving up a hill.**

Forward only. If you reverse while moving up a hill, it may cause the vehicle to tip over.

**NO!**

**Do not attempt to drive across a sloping surface greater than 3°**

Driving across a slope greater than 3° is very dangerous and may cause the vehicle to tip over.

**NO!**

**Do not drive over soft, uneven or unprotected surfaces such as grass, gravel and decks.**


# User Manual

**NO!**

**Do not get on and off on a hill.**


Always stop on the level surface to get in and get out of the vehicle.

**YES!**

**Use low speed while driving down hill.**

When braking while moving down hill, the wheelchair will take longer to come to a complete stop.


**YES!**

Always climb or descend gradients perpendicular to the slope or ramp.

## **Familiarize Yourself With Your Power Wheelchair**

- **Feature Diagram**

### **Kompas**

- |                |                |
|----------------|----------------|
| 1. Handle      | 7. Caster      |
| 2. Backrest    | 8. Footrest    |
| 3. Armrest     | 9. Seat        |
| 4. Fender      | 10. Controller |
| 5. Rear wheel  | 11. Joystick   |
| 6. Drive motor |                |

# User Manual

## ■ Specification

Model No.	Kompas
Seat width	46cm/18"
Seat depth	42cm/16.5"
Seat Height	52cm/20.5"
Overall length	97cm/38"
Overall width (unfolded)	59cm/23"
Overall height	95-100cm/37.5"-39"
Total weight (without batteries)	20.9kg/46lbs
Folded size (L*W*H)	59cm*37cm*81cm/23"*15"*32"
Maximum speed up to	6kph/4mph
Range up to	11km/7.1mi or 18km/11.6mi
Ground clearance	3cm/1.2"
Gradient	6°
Batteries	Std: 24V/11.4AH x 1pc ; Opt: 24V/17.1AH x 1pc
Battery weight (11.4AH/17.1AH)	2.1kg/4.6lbs or 3.0kg/6.6lbs
Motor	DC24V,70W
Controller	PG nVR2
Caster	7"(150*30) Foam filled tire
Drive wheel	8"(200*50) Pneumatic
Brake	Electromagnetic brakes
Capacity	115kg/253lbs

**NOTE:** Recommended temperature for storage and transportation without battery  
-20 ~60 ( -4 °F ~140 °F ).

## ■ Assembly

Figure 1-1

- 1). Ensure that you receive the battery box, and controller.
- 2). Open the power wheelchair.


Figure 1-2

- 3). Pull out the spring pin button. Ensure the pin out from the holder.
- 4). Insert the battery box into the holder.
- 5). Loosen the pin button. Lock the battery.

Figure 1-3

- 6). Insert the controller into the fixing clamp and adjust it in a convenient position.
- 7). Screw the knob to lock the controller.
- 8). The connection wire of the controller is fixed in the way shown in (1-3) : the front end of the connection wire is stuck into the wire buckle, the middle of the connection wire is fixed at the end of the armrest pipe with a wire harness, and the end of the connection wire is inserted into the back cushion.
- 9). Plug in the controller connector to the lower connector, as shown in (1-3), and hide it in the back cushion.

# User Manual


Figure 1-4

10).Please make sure the hook fastened before riding,as shown in (1-4).

**MAKE SURE THE HOOK AND THE PLATE MUST BE IN PERFECT CONTACT.**

## Armrest Width Adjustment


Figure 2-1

- 1). Unlock (A) to lift the armrest.
- 2). Remove the 2 screws (C) to the other two holes.
- 3). Unscrew (B) and move the washer (D) from inside of the frame to outside and put it between armrest and the frame. Screw into the armrest fixation points (B) through the washer (D).


Figure 2-2

- 4). Place Armrest back to fixation points (A).
- 5). Repeat 1-4 on the other side.

## Armrest Height Adjustment


Figure 3-1

- 1). Remove the screws from the armrest.
- 2). Remove (A) and relocate it between armrest and tube.
- 3). Lock screws on the armrest.
- 4). Repeat 1-3 on the other side.

Figure 3-2

## Armrest Position Adjustment


Figure 4-1

- 1). Move one screw from 2<sup>nd</sup> to the 1<sup>st</sup> hole.
- 2). Repeat Step 1 on the other side.

# User Manual

## Detachable Armrest


Figure 5

- 1). Rotate the Lever (A) 180 degree to unlock to detach.
- 2). Repeat 1 on the opposite side.

## Backrest Adjustment



Figure 6

- 1). Remove the screws within the 'Fixing Holes' on both sides of the backrest.
- 2). Push back the backrest by 5 degrees.
- 3). Gently fix screws within the Fixing Holes on one side of the backrest and then proceed to do the same on the other side. Once gently fixed, proceed to tighten the screws on both sides.
- 4). To move the backrest forward, repeat Step 1 to 3 but fix the screws in reverse positions.

Figure 7

## **Operating**

### **■ A. Controls and Indicators**

Figure 8

**P&G nVR2**

It comprises:

1. Battery gauge (fuel meter)
2. Power ON / OFF push-button
3. Joystick
4. Speed adjust button
5. Speed adjust display
6. Horn button

### **■ B. Driving, Steering and Braking**

Press the power ON/OFF push-button. Move the joystick gently forward. The more you push the joystick forward the faster you will travel, but you will not exceed the speed limit imposed by the speed limit control, regardless of the joystick position.

To reverse, pull the joystick to the rear, beyond its centre (neutral) position. The joystick is also your steering lever. Move it to the left and you turn left. Move it to the right and you turn right.

# User Manual

The Power Wheelchair is very responsive and maneuverable. It can turn within its own length. You should avoid making abrupt turns and jerky steering while you are moving at a brisk speed.

If you want to slow down, move the joystick toward its centre position. When you want to stop or have to make an emergency brake application, simply release the joystick. It will automatically return to the centre (neutral) position. At this moment the armature brakes of the electric motors will cut in and bring the Power Wheelchair smoothly to a standstill.

Once the Power Wheelchair is at a standstill, it is automatically locked in its parking position by the electromagnetic safety brake system.

When climbing curbs, even with no curb riders fitted, you will find that your Power Wheelchair will smoothly overcome low to medium high curbs without the need for drastic speed variations, because the micro-computer of the electronic controller automatically compensates for any extra power required.

## Negotiating Gradients

### \*Caution\*

On loose ground (sand, gravel, grass, etc.) The maximum gradient should not exceed 20% (1 in 5).

Always approach an incline directly, not at an angle, and follow it up or down directly, not in a zig-zag movement.

### \*Caution\*

Never attempt to climb or descend an incline which has a slippery or icy surface.

While going uphill or downhill there is no need for drastic corrective joystick movements. The electronic controller ensures that the extra power required for climbing is automatically applied to the motors, so that the selected speed is maintained.


The same applies when you travel downhill. The armature brake and the electromagnetic safety brake give you full control over your powerchair at all times. This unique brake system enables you to descend at a constant slow speed.

You can stop at any time while you are climbing or descending. Simply release the joystick as on the flat. The electromagnetic safety brake will ensure that the powerchair is positively locked and secured in parking position until you are ready to continue your journey.

The stopping distance on slopes might be significantly longer than minimum brake distance from max speed as listed in Tech Data.

## ■ C. Pushing the Powerchair by Hand

Figure 9

If the Power Wheelchair is to be pushed by hand, proceed as follows:

From behind the vehicle, turn the two uncoupling levers (Figure 9) of the electric motors down.

\*Warning\*

1. Never move the uncoupling levers while the vehicle is operating with powered drive.
2. When the electric motors are disengaged the armature brakes and the electromagnetic safety brakes are out of action, so ensure the powerchair is secured with the hand-brakes when it is parked.

When normal powered operation is to be restored, return both uncoupling levers to normal position.

**Important:** The vehicle will not operate with powered drive unless both uncoupling levers are in normal position.

# User Manual

## ■ D. Shutdown

Before you leave your Power Wheelchair, press the ON / OFF push-button. The Warning lamp and fuel meter indicator will extinguish.

Make this shutdown sequence a standard procedure and a habit. It is essential for your personal safety.

## ■ E. Getting Ready for Permanent Use

Now that you have had a good look at your brand new Power Wheelchair and taken it for a test drive in your home and around the house.

It is absolutely necessary to carry out a battery charging procedure before you start normal operation and take your Power Wheelchair on its first long-range outing. Read the next section carefully.

## ■ **Batteries**

The Power Wheelchair has one battery (See General Specification Data) which is accommodated in battery box under the seat of the powerchair. The battery is sealed and maintenance free.

**\*Warning\* Do not attempt to open the batteries.**

Regular monitoring of the battery charge condition and timely recharging are essential to ensure reliability and performance of your Power Wheelchair as well as long battery service life.

The intervals at which the batteries have to be recharged (hours traveled or mileage covered) depend on various factors. It is not possible to specify a generally valid fixed date. Observe the following and you will soon establish a cycle in accordance with your individual requirements and driving routine.

At full battery charge all LED segments of the battery gauge, are lit continuously. With the progressive discharge of the batteries, successive segments will extinguish.

**Note:** fuel meter is only supplied on this model.

If the battery guages show just red and yellow. Depending on age and condition of the batteries, you are then left with a limited driving range before you must recharge the batteries. (approx. 2mi.) To avoid this critical condition, we recommend charging the batteries overnight whenever you have used your Power Wheelchair. It will not only spare your unpleasant situations en route but will also prolong the battery service life.

Should you ignore this warning, while the battery, gauge is on red, and allow the batteries to become almost fully discharged, the electronic controller will automatically shut down the entire electrical system of the Power Wheelchair and bring it to a halt.

This automatic shutdown of the Power Wheelchair, accompanied by fast flashing of the ON / OFF push-button will also occur in the event of any fault in the wiring, the motors, or the controller itself.

# User Manual

**Note:** While climbing a steep gradient, you may find that the battery warning lamp begins to flash. This does not necessarily mean that the batteries are discharged to the critical level. Recheck the warning lamp after you reach level ground.

If during the course of the day, you have used your Power Wheelchair for more than one hour of actual traveling time, it is recommended that you charge the batteries overnighnt. The special-design charger will switch off automatically when the batteries are charged, so it is completely safe and does not use excess electricity regardless of the time the batteries are connected to the charger.

If your Power Wheelchair remains unused for any length of time (approximately four weeks or more), the batteries must be charged. Charge them at least once a month to maintain their charge.

Disposing and recycling of used batteries. Batteries must always be recycled correctly. Do not dispose of them with your trash bin.

For further information, please contact your nearest recycling facility or local authority.

## ■ Battery Charging

The batteries of your Power Wheelchair are different from ordinary car batteries. Only use the charger specifically (HP0060WL2-24V2A) designed for the motive-power type batteries installed.

### Charging procedure


1. Always insert the charger output plug into the charging socket on the controller. Only then connect the charger input cable to an electrical outlet.
2. Follow the instructions on the front panel of the charger for its operation and the meaning of its indicators.
3. The minimum time required for charging varies and depends on battery condition and discharge level. We recommend, however, that the batteries be charged overnight.

**Note:** The specially designed charger ensures that excess power is not consumed regardless of how long it is switched on and connected to the batteries.

4. On completion of the charging cycle, first disconnect the charger input from the electrical outlet and only then disconnect the charger from the controller socket. Do not leave the charger connected to the controller when its input power is disconnected .This will cause the batteries to discharge.

For replacement of the batteries, contact your dealer.

Figure 10

# User Manual

## ■ Lithium Battery Instructions for Use

### A. Charging instructions

- 1.Keep the battery level high. Use it and charge it every day.
- 2.When the two green lights on the controller's power gauge are off, please charge as soon as possible to avoid the damage from the deep discharge of the battery or to avoid the battery exhaustion affects your journey.
- 3.The charging time is controlled at about 8 hours. When the charger indicator turns to green, it means the battery is fully charged. Do NOT over charge for a long time.
- 4.Use the dedicated genuine lithium battery charger only.
- 5.If the battery will be stored for a long time, make its level as 40% (by charging for 2 hours), then disconnect the battery from the device, put it in a cool and dry place. Also charge and discharge every two months to avoid damages.

### B. Warnings

- 1.Do NOT immerse the battery in the water. Place the battery in a dry place when not in use.
- 2.Do NOT put the battery close to a hot and high temperature source, such as fire or heater.
- 3.Do NOT install the positive and negative poles oppositely.
- 4.Do NOT connect the battery to the power socket directly.
- 5.Do NOT throw the battery into fire or heater.
- 6.Do NOT connect the positive and negative poles of the battery with a metal piece.
- 7.Do NOT transport or store the battery together with metal pieces, such as hairpins, necklaces and etc.
- 8.Do NOT beat, throw or step on the battery or etc.
- 9.Do NOT solder the battery and do NOT pierce the battery with nails or other sharp tools.
- 10.Lithium battery can only be discharged between 0°C - 60°C.
- 11.Do NOT charge the battery below 10°C.

## **C. Cautions**

- 1.Do NOT use or place the battery in a high temperature place (sunlight or in a hot car), otherwise it may cause the overheat, catch fire, fail to function, and shorten the life.
- 2.Do NOT use in a strong static electricity or in a strong magnetic field, otherwise it will easily damage the safety device and bring potential dangers.
- 3.If the electrolyte leaked from the battery enters to your eyes accidentally, do NOT rub your eyes. Rinse your eyes with clean water and go to a doctor immediately; otherwise there might be a permanent damage on the eyes.
- 4.If the battery is found as emits odor, heats, changes color or deforms during used, stored or charged. Remove the battery from the device or charger and deactivate it immediately.
- 5.If the electrode is dirty, wipe it with a dry cloth before use, otherwise it may fail to function due to poor contact.
- 6.Abandoned batteries should be covered with insulating paper to prevent fire and explosion.
- 7.To avoid pollution, do NOT discard the battery at will.

# User Manual

## ■ Folding the Power Wheelchair

Figure 11

**Kompas(Unfolded)**

**Kompas(Folded)**

The powerchair can be folded to a neat and compact size within minutes.

1. Loosen the star knob for adjustment of the controller position and slide the controller back as far as it will go. Tighten the star knob.
2. Fold up the footrest.
3. Pull up the lever which under the seat at the front of the chair and press the handle simultaneously. Seat and back will close together to compact stowing size.
4. If the folded Power Wheelchair is stowed or transported in a horizontal position, place it so that the side with controller is uppermost to prevent damage to the controller.

## ■ **Important Safety Information**

No inflation necessary

7"(150\*30) caster 8" drive wheel

Never stand on the footrests while entering or exiting the Power Wheelchair.  
Standing on the footrests could cause the Power Wheelchair to tip over.

The owner should enter or exit the Power Wheelchair only after the power is turned off.

### **Turning Radius**

It is very important, especially when traveling on inclines, to reduce your speed with the control joystick before making a sharp turn. Never try to turn the Power Wheelchair suddenly while traveling at full speed.

Your Power Wheelchair should not be operated on streets except when using the crosswalks. Operate your Power Wheelchair with extreme caution.

### **Anti-Tipper**

Anti-Tipper are standard with the Power Wheelchair. Replacements for these parts are available through your authorized dealer. Anti-Tipper prevent the Power Wheelchair from tipping backwards. We suggest that Anti-Tipper be used on this model Power Wheelchair at all times.

### **Parts**

Special high strength fasteners are used on the powerchair. Replacements for these parts should be ordered through an authorized dealer.

# User Manual

## \*Warnings\*

Never stand on the footrests as it may cause the powerchair to tip over.

Never connect the battery charger to the main supply until after you have connected it to the powerchair.

It is important to develop safe techniques before attempting to negotiate obstacles such as curbs, ramps and inclines. Never try to go up a steep incline without assistance. Avoid uneven and hilly terrain as the Power Wheelchair may tip over. Avoid grades steeper than 10% (1 foot rise in a ten foot distance) without assistance.

Control your speed carefully when traveling downhill. Slow speeds should be used when traveling on downgrades.

Do not use the ON / OFF switch to turn off the power to the Power Wheelchair while it is moving. Using this switch while the Power Wheelchair is still in motion causes it to come to an abrupt stop and may cause injury to the owner or damage to the Power Wheelchair.

Many power wheelchair controllers have been affected by high power radio transmitters, such as CB radios and cellular telephones. If any of these devices are being used near your powerchair, you should use caution in operation. If loss of control occurs, shut off your Power Wheelchair immediately and leave it off until the interfering devices are no longer operating.

The user should be trained by a healthcare professional to assure safe operation of the powerchair during the performance of daily activities.


## Weight Limit

Maximum weight limit should not exceed standards (See General Specification Data).

## Controller

The electronic controller is what you use to operate your Power Wheelchair. It takes the battery voltage and sends it to the appropriate system. The electronic controller enables you to move the Power Wheelchair, as well as monitor battery status, electronic controller functions, and the conditions of your electrical system. The controller is an integral electronic controller. All of the electronics necessary to operate the Power Wheelchair are contained in one module.

Figure 12

The controller consists of:

1. Joystick
2. Keypad
3. Off-board charger / programming socket
4. Motor connector

Typically, the controller is mounted to one of the armrests and is connected to the motors and batteries.

## Joystick

The joystick controls the direction and speed of your Power Wheelchair. When you move the joystick from the neutral position (center) the electromagnetic brakes release and allow your Power Wheelchair to move. The further you push the joystick from its neutral position, the faster the Power Wheelchair will move. When you release the joystick and allow it to return to the neutral position, you engage the electromagnetic brakes. This causes the Power Wheelchair to decelerate and come to a complete stop.

**If your Power Wheelchair begins to move in an unexpected manner, immediately release the joystick. Unless the joystick is damaged, this action should stop your Power Wheelchair.**

# User Manual

## P&G nVR2 Controller

### Keypad

The keypad is located on the front of the joystick.

It contains keys necessary to operate your power chair.

Figure 13

**P&G nVR2**

### 1. On / Off Button

The on / off button turns the controller on and off.

**Unless faced with an emergency situation, do not use the on / off key to stop the chair. This will cause the power chair to stop abruptly.**

**Always turn the power off when you are stationary to prevent unexpected movement.**

### 2. Battery Condition Meter

#### P&G nVR2:

The battery condition meter is located on the front of the joystick. This is a 10-segments illuminated display that indicated that the nVR2 is turned on and also gives the battery status, the nVR2 status, and the electrical system status.

**1 LED**  
The battery needs charging or there is a bad connection to the battery. Check the connections to the battery. If the connections are good, try charging the battery.

**2 LED**  
The left hand motor\* has a bad connection. Check the connections to the left hand motor.

**3 LED**  
The left hand motor\* has a short circuit to a battery connection. Contact your service agent.

**4 LED**  
The right hand motor\* has a bad connection. Check the connections to the right hand motor.

**5 LED**  
The right hand motor\* has a short circuit to a battery connection. Contact your service agent.

**6 LED**  
The wheelchair is being prevented from driving by an external signal. The exact cause will depend on the type of wheelchair you have, one possibility is the battery charger is connected.

**7 LED**  
A joystick fault is indicated. Make sure that the joystick is in the center positon before switching on the control system.

**8 LED**  
A possible control system fault is indicated. Make sure that all connections are secure.

**9 LED**  
The parking brakes have a bad connection. Check the parking brake and motor connections. Make sure the control system connections are secure.

**10 LED**  
An excessive voltage has been applied to the control system. This is usually caused by a poor batteryconnection. Check the battery connections.

**7 LED +S**  
A communication fault is indicated. Make sure that the joystick cable is securely connected and notdamaged.

**8 LED +A**  
An Actuator trip is indicated. If more than one actuator is fitted, check which actuator is not workingcorrectly. Check the actuator wiring.

# User Manual

### 3. Speed / Profile Buttons

There are two buttons that control either the speed or the profile. This depends on how your controller was programmed. Press the speed / profile increase button to increase the speed or change the profile. Press the speed / profile decrease button to decrease the speed or change the profile.

The speed / profile setting is displayed on the maximum speed / profile indicator. If your powerchair was programmed with a drive profile, contact your authorized dealer to provide more information.

**NOTE:** We recommend that the first few times you operate your powerchair, you set the speed to the slowest setting until you become familiar with your new powerchair.

### 4. Horn Button

The horn button activated the horn.

#### Off - Board Charger / Programming Socket

You may use an off - board charger to charge the powerchair batteries through the 3 - pin socket located in the front of the controller. If you use an off - board charger, the charger current should not exceed 12 amps. contact your dealer for more information.

**Note:** The socket may also be used for programming the controller. Contact your dealer for more information.

**NOTE:** If any of the above LEDs flash rapidly, it means either a total short circuit, a single open circuit, or a total open circuit in the associated indicator has been detected.

**CAUTION!** Failure to properly align the connectors can result in damage to the controller, the charger, and the connectors.

#### Battery Connector

This connects the controller to the powerchair's battery box.

#### Motor Connector

This connects the controller to the powerchair's motors and brakes.

## **Thermal Rollback**

The controller is equipped with a thermal rollback circuit. The circuit monitors the temperature of the controller, which roughly translates to motor temperature. In the event that the controller becomes excessively hot (above 140 deg. F or 60 deg. C) motor current (amperage) is reduced. For every degree above 140 deg. F, the motor current limit is reduced by 40 amps until the controller reaches 158 deg. F or 70 deg. C, at which time the current output is reduced to zero. This reduces your chair's "power", which also could reduce your chair's speed, and allows the electrical components and motors to cool down. When the temperature returns to a safe level, your powerchair resumes normal operation.

## **VR2 Error Codes**

The VR2 controller is designed with the user's safety as the prime consideration. It incorporates many sophisticated self-test features which search for potential problems at a rate of 100 times per second.

If the VR2 detects a problem either in its own circuits or in the powerchair's electrical system, it may stop the Power Wheelchair, depending on the severity of the problem. The VR2 is designed to maximize the user's safety under all normal conditions. The table below identifies the individual error codes. Error codes are displayed as a rapid flashing of lights. If you get one of these error codes, please contact your dealer.

# User Manual

Flashing Lights	Diagnosis / Solution
1	The battery needs charging, or there is a bad connection to the battery. Check the connections to the battery. If the connections are good, try changing the battery.
2	The left motor has a bad connection. Check the motor connection.
3	The left motor has a short circuit to a battery connection. Contact your dealer.
4	The right motor has a bad connection. Check the motor connection.
5	The right motor has a short circuit to a battery connection. Contact your dealer.
6	The power chair is being inhibited by the battery charger. Unplug the battery charger.
7	A joystick fault is indicated. Make sure that the joystick is in the neutral (center) position before turning on the controller.
8	A controller fault is indicated, make sure that all connections are secure.
9	An excessive voltage has been applied to the controller system. This is usually caused by a poor battery connection. Check the battery connections.

## Battery Boxes

On the top right side of the rear battery box, you will see the reset button of a circuit breaker. This button must always be pressed in, otherwise the powerchair will not operate. \*Refer to the Troubleshooting Guide\* Also on the battery boxes are cables with connectors for connecting the batteries together and to connect batteries to the controller box.

## ■ Maintenance

There are no specific maintenance requirements for the joystick and the controller, but the service life of the equipment will be extended if you observe the following:

The joystick box and the controller are splash proof but should be protected from excessive spillage of liquid. Do not expose to direct heat (from an electric heater, for example). Check that the rubber boot on the joystick and switch cover are in good condition and show no holes or cracks.

### **Wheels - Removal / Installation**

#### **A. Rear Wheel**

1. Ensure that the freewheel lever of the motor drive is engaged. (Figure 4)
2. Support or jack up the powerchair under the frame side-member.
3. Unscrew the hexagon nut and remove, together with the washer.
4. Withdraw the rear wheel from the drive axle.

Installation is a reversal of the removal procedure. Check that the threads of the hexagon nut are undamaged and clean. If necessary, wipe with a clean dry cloth. Do not use any sharp tools to clean the threads. Do not forget to slide the washer over the nut before it is screwed in. Tighten the nut firmly.

#### **B. Front Wheel**

1. Support or jack up the Power Wheelchair under the frame side-member.
2. Using the wrench provided, remove the axle bolt.
3. Slide the front wheel out of the fork.

Installation is a reversal of the removal procedure.

#### **C. General Cleaning**

The shrouds and seat of the wheelchair should be cleaned using a damp cloth and gentle detergent. Do not use abrasive cleaners or solvents, as this will damage the plastic components. Do not spray water (hose or pressure washer) onto the wheelchair as this may cause damage to the electronics.

#### **D. Long-term storage**

The packaging of the wheelchair should be stored for a further storage or transport that might become necessary.

# User Manual

## ■ Troubleshooting Guide

Nature of Fault	Problem Cause	Action
1.Powerchair will not move, power on lamp doesn't light up.	a) Master plug not inserted b) ON / OFF push-button in off position c) Batteries discharged to low level d) Loose plug connection at controller or battery	a) Insert master plug b) Press in c) Re-charge the batteries d) Press plugs firmly into sockets
2.Powerchair will not move, power on lamp lights up.	a) Drive motor not engaged b) Fault in controller electronics	a) Reset uncoupling lever(s) b) Consult dealer

### \*Important\*

In all cases of an electrical fault, check that the circuit breaker (located on the rear battery box) is pressed in. If it has been tripped it must be pressed. If the reset button comes out again, even after the suspected cause of the fault has been corrected, consult your dealer.

## IEC SYMBOLS

	<b>Caution, attention or consult accompanying documents.</b>
	<b>Alternating Current</b>
	<b>Type BF Equipment</b>
	<b>Double Insulation</b>
	<b>No Smoking or Naked Flames</b>
	<b>Not intended to use as a seat in a motor vehicle</b>

Degree of protection against ingress of water is rated as IPx0.

Not intended to use as a seat in a motor vehicle.

# User Manual

Disclosure information(ISO)					
Standard reference	min.	max.	Standard reference	min.	max.
Overall length with legrest	—	970mm	Seat plane angle	—	5°
Overall width	—	590mm	Effective seat depth	—	420mm
Folded length	—	590mm	Effective seat width	—	460mm
Folded width	—	370mm	Seat surface height at front edge	—	520mm
Folded height	—	810mm	Backrest angle	—	12°
Total mass	23kg	23.9kg	Backrest height	—	330mm
Mass of the heaviest part	2.1kg	3.0kg	Footrest to seat distance	420mm	480mm
Static stability downhill	—	6°	Leg to seat surface angle	—	95°
Static stability uphill	—	6°	Armrest to seat distance	—	220mm
Static stability sideways	—	3°	Front location of armrest structure	—	—
Energy consumption	11km	18km	Hand rim diameter	—	—
Dynamic stability uphill	—	6°	Horizontal location of axle	76mm	100mm
Obstacle climbing	—	50mm	Minimum turning radius	—	825mm
Maximum speed forward	—	6km/h			
Minimum braking distance from max speed	—	1000mm			

## ■ **Limited Warranty**

The terms and conditions of the warranty are part of the general terms and conditions particular to the individual countries in which this product is sold.

# User Manual

We wish you a safe and comfortable riding experience!

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Invacare AG
Invacare Poirier SAS
Merits Health Products Co., Ltd.
MedNet EC-REP GmbH
Invacare GmbH