
Flo-tech® Xtra Accessories

en	User Manual	2
de	Gebrauchsanweisung	4
da	Brugsanvisning	6
es	Manual del usuario	9
fi	Käyttöohje	11
fr	Manuel d'utilisation	13
it	Manuale d'uso	16
nl	Gebruiksaanwijzing	18
no	Bruksanvisning	20
pt	Manual de utilização	23
sv	Bruksanvisning	25

This manual MUST be given to the user of the product.
BEFORE using this product, this manual MUST be read and saved for future reference.

Yes, you can.®

1 General

1.1 Introduction

This user manual contains important information about the handling of the product. To ensure safety when using the product, read the user manual carefully and follow the safety instructions.

Invacare reserves the right to alter product specifications without further notice.

Before reading this document, make sure you have the latest version. You find the latest version as a PDF on the Invacare website.

If you find that the font size in the printed document is difficult to read, you can download the PDF version from the website. The PDF can then be scaled on screen to a font size that is more comfortable for you.

In case of a serious incident with the product, you should inform the manufacturer and the competent authority in your country.

1.2 Compliance

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485.

This product features the CE mark, in compliance with the Medical Device Regulation 2017/745 Class I.

This product features the UKCA mark, in compliance with Part II UK MDR 2002 (as amended) Class I.

We are continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

We only use REACH compliant materials and components.

1.3 Warranty Information

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business in the respective countries.

Warranty claims can only be made through the provider from whom the product was obtained.


1.4 Limitation of Liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorised modifications and/or use of unsuitable spare parts

2 Product Overview

2.1 Product Description

(A)	Pelvic obliquity pad
(B)	Clinical fluid sac
(C)	Fluid pack
(D)	Side wedges
(E)	Lateral wedges
(F)	Front wedge
(G)	Deep pommel
(H)	Standard pommel
(I)	Thigh aligner
(J)	Hip aligner

2.2 Intended Use

These accessories are intended to be used in conjunction with a compatible cushion to meet the specific positioning needs of the end user. Cannot be used on their own.

Intended users

Flo-tech accessories can be used for users who need pelvic stability and are at very high risk of developing pressure ulcers.

Indications / contra-indications

There are no contra-indications known.


3 Accessories / Options

3.1 Installing accessories


All the parts are installed using the hook and loop strips.

ALWAYS remove the clinical fluid sac (B) before fitting the accessories.


Installing Fluid pack

1. Remove the outer cover, inner liner and clinical fluid sac ⑧.
2. Align the hook strip under the fluid pack ③ with the loop on the base of the cushion.
3. Press down to secure the fluid pack to the base of the cushion.
4. Re-install the clinical fluid sac, inner liner and outer cover.

Installing Abductor Support — Standard pommel or Deep pommel


1. Remove the outer cover, inner liner and clinical fluid sac.
2. Align the hook strip under the abductor support ⑥ or ⑦ with the loop on the base of the cushion.
3. Press down to secure the abductor support to the base of the cushion.
4. Re-install the clinical fluid sac, inner liner and outer cover.

Installing Hip and Thigh Aligners

1. Remove the outer cover, inner liner and clinical fluid sac.
2. Position the thigh aligner ① and/or the hip aligner ① on the base of the cushion aligning the hook and loop strips and press down firmly.
3. Re-install the clinical fluid sac, inner liner and outer cover.

Installing Pelvic Obliquity Pad

1. Remove the outer cover, inner liner and clinical fluid sac.
2. Position the pelvic obliquity pad ① on the base of the cushion aligning the hook and loop strips and press down firmly.
3. Re-install the clinical fluid sac, inner liner and outer cover.

Installing Wedges

1. Remove the outer cover, inner liner and clinical fluid sac.
2. Position the wedges in the desired position on the underside of the base of the cushion as shown and press down firmly.
 - Side wedge ④
 - Lateral wedge ⑤
 - Front wedge ⑥
3. Re-install the clinical fluid sac, inner liner and outer cover.

 Foam wedges can be trimmed to size as needed.

Trimming the Xtra parts and base

The unique components of the Xtra Modular System can be trimmed to provide the ideal seating position for the end user. The front, lateral and side wedges require trimming to fit the base. These parts are intended to indicate the size of the cushion and act as a marker to help trim the part accurately. Use a precision knife or scissors to trim the accessories (please see diagrams for examples). The base of the cushion can also be trimmed to accommodate leg length discrepancies, additional hamstring space for knee flexion contractures or amputations.

4 Maintenance

4.1 General Maintenance Information

WARNING! Risk of Injury

- Maintenance should be done at least once a year.
- Do not use a defective product.
- Do not attempt repairs yourself. For repairs contact your provider.

4.2 Inspection

Check the accessories for strike-through (this may include fluid ingress, stains, rips or damage) after the use of each patient and on a minimum monthly basis (depending on which occurs first) by a suitably qualified and competent person.

Check accessories

1. Unzip the cover completely.
2. Check for any staining on the interior foam.
3. Replace any damaged cover or stained foam and dispose of as per local authority procedure.

4.3 Cleaning and Disinfection

General Safety Information

CAUTION! Risk of Contamination

- Take precautions for yourself and use appropriate protective equipment.

NOTICE!

- Wrong fluids or methods can harm or damage the product.
- All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.
 - Never use corrosive fluids (alkalines, acid etc.), abrasive cleaning agents or solvents (cellulose thinner, acetone etc.). We recommend a mild detergent.
 - Always make sure that the product is completely dried before taking it into use again.

For cleaning and disinfection in clinical or long-term care environments, follow your in-house procedures.

Cleaning Intervals

NOTICE!

Regular cleaning and disinfection enhance smooth operation, increases the service life and prevents contamination.

Clean and disinfect the product:

- regularly while in use,
- before and after any service procedure,
- when it has been in contact with any body fluids,
- before using it for a new user.

5 Technical Data

5.1 Materials

Components	Type of material
Fluid sac Fluid pack	Polyether based thermoplastic polyurethane film with silicone gel
Pelvic obliquity pad Side wedges Lateral wedges Front wedge Deep pommel Standard pommel Thigh aligner Hip aligner	High resilient polyurethane foam

All product components do not contain natural rubber latex.

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Invacare Limited
Invacare UK Operations Limited