Invacare® Dacapo®

Dacapo® Basic, Dacapo® Basic Low, Dacapo® 2Sided,
Dacapo® Premier, Dacapo® Kid, Dacapo® Excel,
Dacapo® Ampla, Dacapo® Top, Dacapo® Top Pro

en Mattress range

User Manual 2

da Madrasserie

Brugsanvisning 7

de Matratzen

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es Gama de colchones

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fi Patjavalikoima

Käyttöohje. 27

fr Gamme de matelas

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it Gamma di materassi

Manuale d'uso. 39

nl Matrassen

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no Madrassområde

Bruksanvisning 52

pt Gama de colchões

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sv Madrassortiment

Bruksanvisning 64

This manual MUST be given to the user of the product.
BEFORE using this product, this manual MUST be read and
saved for future reference.

Yes, you can.®

1 General

1.1 Introduction

This user manual contains important information about the handling of the product. To ensure safety when using the product, read the user manual carefully and follow the safety instructions.

Only use this product if you have read and understood this manual. Seek additional advice from a healthcare professional who is familiar with your medical condition and clarify any questions regarding the correct use and necessary adjustment with the healthcare professional.

Note that there may be sections in this document, which are not relevant to your product, since this document applies to all available models (on the date of printing). If not otherwise stated, each section in this document refers to all models of the product.

The models and configurations available in your country can be found in the country-specific sales documents.

Invacare reserves the right to alter product specifications without further notice.

Before reading this document, make sure you have the latest version. You find the latest version as a PDF on the Invacare website.

Previous product versions may not be described in this manual's current revision. If you require assistance, please contact Invacare.

If you find that the font size in the printed document is difficult to read, you can download the PDF version from the website. The PDF can then be scaled on screen to a font size that is more comfortable for you.

For more information about the product, for example product safety notices and product recalls, contact your Invacare distributor. See addresses at the end of this document.

In case of a serious incident with the product, you should inform the manufacturer and the competent authority in your country.

1.2 Symbols in this Document

Symbols and signal words are used in this document and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

WARNING!

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

CAUTION!

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

NOTICE!

Indicates a hazardous situation that could result in damage to property if it is not avoided.

Tips and Recommendations

Gives useful tips, recommendations, and information for efficient, trouble-free use.

Other Symbols

(Not applicable for all manuals)

Triman

Indicates recycling and sorting rules (only relevant for France).

UKRP

UK Responsible Person. Indicates if a product is not manufactured in the UK.

1.3 Compliance

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485.

This product features the CE mark, in compliance with the Medical Device Regulation 2017/745 Class I.

This product features the UKCA mark, in compliance with Part II UK MDR 2002 (as amended) Class I.

We are continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

We only use REACH compliant materials and components.

1.4 Warranty Information

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business in the respective countries.

Warranty claims can only be made through the provider from whom the product was obtained.

1.5 Limitation of Liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorised modifications and/or use of unsuitable spare parts

1.6 Service Life

The expected service life of this product is five years when used daily and in accordance with the safety instructions, maintenance intervals and correct use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

2 Safety

2.1 Safety Information

Education, clinical judgement and action-based planning based on vulnerability are fundamental factors in prevention of pressure injuries.

A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

WARNING!

Risk of Serious Injury or Damage

Improper use of this product may cause injury or damage.

- If you are unable to understand the warnings, cautions or instructions, contact a health care professional or provider before attempting to use this equipment.
- Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manual, service manual or instruction sheet supplied with this product or optional equipment.

WARNING!

Risk of Fire or Explosion!

A cigarette can burn a hole in the bed surface and cause damage to the device. Also, patient clothing, bed sheets, etc. may be combustible and cause a fire. Failure to observe this warning can result in a severe fire, property damage and cause physical injury or death.

- Take special care in oxygen rich environments.
- Do not smoke.

WARNING!

Risk of Developing Pressure Injuries

The device delivers effective pressure redistribution, when the support surface is covered with a cotton, cotton combination or linen bed sheet, and any one of these is the only item deployed between the support surface and the user.

- Bed sheets must be loosely fitted, with creases smoothed out.
- Ensure that the support surface in contact with the user is kept free from crumbs and other food debris, and that drip cables, stents, and other foreign objects do not become entrapped between the support surface and the user.
- Heated over blankets must only be used in consultation with a qualified healthcare professional, as an increase in temperature can increase the risk of developing pressure injuries.

CAUTION!

Risk of Injury

Non-original or wrong parts may affect the function and safety of the product.

- Only use original parts for the product in use.
- Due to regional differences, refer to your local Invacare catalogue or website for available options or contact your Invacare distributor. See addresses at the end of this document.

CAUTION!

Risk of Injury or Damage to Property

- Mattress must be used with the cover at all times.

2.2 Safety Information for Transport

- Take care when handling the product to ensure no damage to the cover.
- It is recommended that two people lift/carry the product.
- Avoid contact with jewellery, nails, abrasive surfaces, etc.
- Do not drag the product.
- Avoid contact with wall, door frames, door catches or locks, etc.
- Do not transport in roll cages unless completely protected from the sharp edges of the cage.

3 Product Overview

3.1 Intended Use

This mattress is intended to deliver effective pressure reduction and postural support to users. It is intended to be used in conjunction with an appropriately sized bed frame as part of an overall pressure injury prevention program of care.

Indications

Suitable for users with a low risk of developing pressure injuries and users who require a comfortable mattress. **Dacapo 2Sided, Premier, Ampla** and **Excel** are pressure injury-preventive and can be used in cases of medium to high risk of developing pressure injuries and for the care of pressure injuries at stage 1 or 2.


3.2 Symbols on the Product

	European Conformity		Medical device
	Manufacturer		UK Conformity Assessed
	European representative		Date of manufacture
	Importer		LOT number
	Do not pierce or cut		Unique device identifier
	Maximum user weight. See Technical data.		Read user manual
	Do not put near flame		Dry clean, PCE only
	Do not iron		Do not bleach

	Machine wash temperature. Max. 95 °C		Tumble dry low heat
	Max. 40 °C (Cotton cover)		
	Line dry		Orientation symbol

3.3 Scope of Delivery

The following components are included within the scope of delivery:

A	Mattress with cover
B	User Manual

3.4 Mattresses

This range of mattresses offers pressure reducing qualities, specific to the individual's needs. By promoting even weight distribution, they reduce tissue pressure at potentially vulnerable areas.

Dacapo Basic

Its high density foam provides a firm supportive base. Fully satisfies basic comfort and hygiene requirements for a sleeping surface.

Dacapo Basic Low

Its high density foam provides a firm supportive base.

Dacapo 2Sided

Multi-layered support system, with a high resilience foam on one side and viscoelastic foam on the other side. The high resilience side has been developed with zoned castellation. The viscoelastic side adapts itself both to the user's body dynamics and temperature.

Dacapo Premier

Features a unique contoured insert pad which provides patient comfort and optimises pressure reduction. The insert is supported by a U-core foam that facilitates patient transfer by providing effective side-wall support.

Dacapo Kid

Constructed from a single layer of high resilience foam. Its unique square castellation, cut into the mattress, minimises surface pressure, which allows pressure redistribution.

Dacapo Excel

Multi-layered support system of high resilience foam. The top layer's zoned castellation minimises surface pressure, which allows for optimal pressure redistribution.

Dacapo Ampla

Single-layered support system which offers firm responsive support, pressure redistribution and comfort. It is designed for the plus size user.

Dacapo Top

Overlay designed to be placed on top of basic mattresses. It is composed of viscoelastic foam which adapts to the user's body dynamics and temperature.

Dacapo Top Pro

Overlay designed to be placed on top of basic mattresses. It is composed of a castellated high resilience foam which provides comfortable support.

3.5 Mattress Covers

The Dacapo range has three cover options available:

Wipedown Cover

The wipedown cover is composed of polyurethane coated polyester. The cover is two-way stretch and gentle on the skin. The fabric is vapour-permeable and water-resistant.

This cover is allergy-tested.

Evacuation Cover

The evacuation cover is composed of two-way stretch, vapour-permeable, water-resistant, polyurethane coated polyester. These covers are equipped with handles on all four sides. On the two long sides there are straps, which in case of an emergency evacuation can be wrapped around the user.

This cover is allergy-tested.

Cotton Cover

The cotton cover is composed of standard cotton, washable at 40 °C.

4 Usage

4.1 Safety Information

WARNING!

Insufficient repositioning of the patient can lead to tissue compression and potential injury formation. To relieve pressure, it is very important for the patients to reposition themselves, or to be repositioned, on a regular basis.

- Before using the product, always consult a qualified healthcare professional for clinical judgement.
- Monitor the patient frequently.

**CAUTION!**

- Make sure that the printed side of the mattress cover always faces upwards.
- Ensure that the distance between the surface of the mattress and the top of the side rail is at least 220 mm.

**CAUTION!****Risk of damage to the mattress**

If holes are present in the mattress cover, there is a risk that liquids may ingress and contamination may occur.

- Ensure that the mattress is not jammed or damaged by sharp edges.
- Do not place hypodermic needles, venflons, scalpels or other similarly sharp objects on or under the mattress.
- Do not use electrically heated blankets directly on or under the mattress.
- Ensure that all venflons are taped down correctly with no sharp edges exposed.
- When using bridging boards or other patient transfer aids, check for any sharp edges or burrs before use. When using the mattress on a profiling bed ensure that the knee break is used before the backrest.
- Attach medical equipment including infusion pumps and monitors to appropriate bed accessories.
- Avoid cigarette burns and pet claws puncturing the mattress cover.

4.2 Preparing the Product for Use

**CAUTION!****Risk of Injury or Damage to Property**

- Before use, check the product for any defect.
- Do not use a defective product.

1. Remove all packaging before use.
2. Place the mattress directly on the mattress support of the bed or place the mattress overlay on the existing mattress.

Compress rolled products must be left to recover before use.

Follow the instructions on the packaging inlay.

Products that are delivered flat-packed can be used immediately.

4.3 Using the Handles

Handles on the Wipedown Cover

**CAUTION!****Risk of Injury**

The handles on the cover are meant to transport the mattress only.

- Do not use the handles to move patients on the mattress.

Handles on the Evacuation Cover

The evacuation cover is equipped with handles on all four sides. On the two long sides are straps, which can be wrapped around the patient.

When the straps are not in use, they can be tucked away under the hook-and-loop pads. The handles can be used for lifting/dragging the mattress.

1. In case of an emergency evacuation, wrap the straps on the long sides around the patient.
2. Make sure that the arms and legs are within the straps, and the straps are tight so that the mattress will bend up around the patient.

5 Maintenance

5.1 Inspection

Check product (foam and cover) for strike-through (this may include fluid ingress, stains, rips or damage) after the use of each patient and on a minimum monthly basis (depending on which occurs first) by a suitably qualified and competent person.

Check Overlays

1. Unzip the cover completely.
2. Check for any staining on the white underside of the cover.
3. Check for any staining on the interior foam.
4. Replace any damaged cover or stained foam and dispose of as per local authority procedure.

5.2 Cleaning and Disinfection

General Safety Information

**CAUTION!****Risk of Contamination**

- Take precautions for yourself and use appropriate protective equipment.

**NOTICE!**

Wrong fluids or methods can harm or damage the product.

- All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.
- Never use corrosive fluids (alkalines, acid etc.) or abrasive cleaning agents. We recommend an ordinary household cleaning agent such as dishwashing liquid, if not specified otherwise in the cleaning instructions.
- Never use a solvent (cellulose thinner, acetone etc.) that changes the structure of the plastic or dissolves the attached labels.
- Always make sure that the product is completely dried before taking into use again.

For cleaning and disinfection in clinical or long-term care environments, follow your in-house procedures.

Cleaning Intervals

- ! NOTICE!**
Regular cleaning and disinfection enhance smooth operation, increases the service life and prevents contamination.
Clean and disinfect the product
- Regularly while it is in use,
 - Before and after any service procedure,
 - When it has been in contact with any body fluids,
 - Before using it for a new user.

- ! NOTICE!**
- Keep a cleaning record as part of cleaning the system.

Cleaning Instructions

Cleaning Covers

(Removal of contaminants such as dust and organic matter)

1. Remove cover for laundering.
2. Launder the cover with maximum temperature as stated on product label, using a diluted detergent solution (Instructions on label).

- ! NOTICE!**
- Washing at higher temperatures will cause shrinkage.

Drying Covers

1. Hang cover from a line or bar and drip dry in a clean indoor environment or;
2. Tumble dry on a low temperature.

- ! NOTICE!**
- Tumble dry setting must not exceed 40 °C.
 - Do not tumble dry for longer than 10 minutes.
 - Dry thoroughly before use.

Disinfection Instructions

- ! NOTICE!**
- Only use disinfectants and methods approved by your local infection control institution and follow your local infection control policy.
 - Follow your local decontamination protocols.

Disinfecting Covers

(Reducing the number of microorganisms)

- ! NOTICE!**
Failure in disinfection process may result in the accumulation of reagent that could damage the polyurethane coating, react with the bed frame, or negate the biocompatibility results.
- Ensure that all cleaning agents, and disinfectants, are thoroughly rinsed off.
 - Dry thoroughly before use.

1. Wipe down the cover with a suitable detergent.
2. Rinse the cover thoroughly with clean water using a single use nonabrasive cloth.
3. Dry the cover thoroughly.

- ! NOTICE!**
Heavy Soilage
Where the mattress is badly soiled, we recommend cleaning with a dilute cleaning solution at maximum temperature in the washing machine (see product label).
- Clean up all spillages of bodily fluids e.g., blood, urine, faeces, sputum, wound exudate and all other bodily secretions as soon as possible using a suitable detergent.
 - Large spillages of blood should be absorbed and removed with paper towels first.

- WARNING!**
- In event of contamination, contact your hygiene specialist.
 - Remove contaminated foams from use.

Autoclaving Covers

Autoclave the cover at 110 °C.

- ! NOTICE!**
- Make sure the cover is outstretched (not folded).
 - Do not place covers on top of each other.

5.3 Replacing Cover

1. Unzip cover and remove it carefully from the foam core.
2. Place new cover onto the foam core.
3. Close zipper.

- ! NOTICE!**
- Ensure that the corners of the foam core are positioned correctly into the corners of the cover.
 - Ensure that the castellated foam is facing upwards when inserted into its cover.

- WARNING!**
- In event of contamination, contact your hygiene specialist.
 - Remove contaminated foams from use.

6 After Use

6.1 Storage

- ! NOTICE!**
- Store device in a dry environment.
 - Store device within a protective cover.
 - Carefully roll up the device and store it on a clean, dry (off-floor) surface with no sharp edges to prevent possible damage.
 - Never store other items on top of a device.
 - Do not store device next to radiators or other heating devices.
 - Protect device from direct sunlight.

6.2 Reconditioning

This product is suitable for reuse. To recondition the product for a new user, carry out the following actions:

- Inspection according to service plan
- Cleaning and disinfection

For detailed information, see *5 Maintenance*, 5.

Make sure that the user manual is handed over with the product.

If any damage or malfunction is detected, do not reuse the product.

6.3 Disposal

Be environmentally responsible and recycle this product through your recycling facility at its end of life.

Disassemble the product and its components, so the different materials can be separated and recycled individually.

The disposal and recycling of used products and packaging must comply with the laws and regulations for waste handling in each country. Contact your local waste management company for information.

7 Technical Data

7.1 General Data

Dimensions

Product	Dimensions [mm]		
	Width	Length	Height
Dacapo Basic	800, 850, 900	2000	120
Dacapo Basic Low	830, 900	2000	60
Dacapo 2Sided	800, 850	2000, 2050	150
Dacapo Premier	800 - 1200	2000, 2100	150
Dacapo Kid	600 - 800	1358 - 1800	100
Dacapo Excel	800 - 1200	2000, 2100, 2200	120, 140
Dacapo Ampla	850 - 1200	2000	150
Dacapo Top	750 - 1200	1800 - 2200	70
Dacapo Top Pro	800 - 1200	2000, 2150, 2200	70

Weights

Product	Maximum user Weight [kg]	Weight of Product [kg] ¹
Dacapo Basic	130	8.5 - 12
Dacapo Basic Low	130	3.5 - 5
Dacapo 2Sided	130	11.5 - 17
Dacapo Premier	180	14
Dacapo Kid	70	5 - 6.5
Dacapo Excel	130, 150	10 - 15
Dacapo Ampla	350	10 - 14
Dacapo Top	130	5 - 7.5
Dacapo Top Pro	150	4.5 - 5

¹ Weight can vary depending on product size. Average weight used as an indication.

7.2 Materials

Foam	Polyurethane Foam
Wipedown cover / Evacuation cover	Polyurethane coated polyester knitted fabric
Cotton cover	100 % cotton fabric

All product components do not contain natural rubber latex.

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Invacare AB
Invacare UK Operations Limited