


## Data sheet of the e-fix E35/E36 lithium ion battery in terms of transportation



Artikel	1565639	/	1566153
Rated cell capacity:	2,0 Ah	/	2,5 Ah
Rated battery capacity:	6,0 Ah	/	7,5 Ah
Number of cells:	30	/	30
Watt-hour rating:	216 Wh	/	270 Wh

Both, cells and the battery pack are tested according the UN Manual "Test and Criteria" part III, 38.3.

We hereby certify that the lithium ion battery pack of the e-fix E35/E36 meets the requirements of the transportation as **dangerous goods class 9**

### UN 3480

under the regulations of **ADR, RID, IMDG and IATA**.

Condition precedent to meet the requirements is the transportation in the tested dangerous goods carton (see UN number on the carton) provided by Ulrich Alber.

Due to frequent amendments of the respective regulations this document is valid until 31<sup>st</sup> December 2015.

Albstadt, 11<sup>th</sup> March 2015

Alber GmbH



Ralf Ledda  
Managing Director