

ZERTIFIKAT ◆ CERTIFICATE ◆ 認証証書 ◆ CERTIFICADO ◆ CERTIFICAT

Certificate

No.: Z1 02 09 42940 006



Invacare EC-Hong A/S

Ostergade 3
4270 Hong
Denmark

with production facility(ies)
42940

is authorized to label the following products with the
TÜV Mark P.
Please see notes overleaf.

Product: Pflegebetten
Patient Care Beds

Model: Patient Care Bed
Type: SB 400

Parameters:

mechanical data	
bed:	1000 x 2130 mm
height of the bed side min./max.:	700/1260 mm
height of the mattress support	
platform min./max.:	330/900 mm
safe working load max.:	180 kg
patient weight max.:	145 kg

electrical data	
connection voltage:	230 V, 50 Hz
power consumption max.:	1.5 A

controller: co. Linak, 24 Volt, max. 120 VA,
IP 66, type CB09L0-3t-24
adjustment of the back section: co. Linak, 24 Volt,
max. 3,8 A, IP 66, type LA31.2J-130-24-301, 4000 N
level adjustment: co. Linak, 24 Volt, max. 3,8 A,
IP 66, type LA-31.40BM-250-24-001

The compliance of the above-mentioned product with the following essential requirements
was tested on a voluntary basis:

EN 1970:2000

Report No: 04610045159, E - (Z1 99 02 10601 015)

Released with the above certificate number by the
certification body of TÜV PRODUCT SERVICE GMBH.

Department: HANMED5 / zr
Date: 24-Sep-2002

