



Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that our products:

e-fix (E20, E25, E26)^{P(small),M(extended)}
viamobil (V15), viamobil eco (V14)^P
viaplus (V12)^P
scala (S30, S31, S35, S36, S38, S39)^P
scala eco S34^{*P}
Adventure (A10)^{P,M}

including the non-spillable valve regulated lead acid batteries (VRLA) fulfilling special provision A67 (therefore not restricted according to UN 2800) are appropriate to be carried by passengers according

IATA chapter 2.3.2.2

and therefore may be carried as checked baggage or in the cargo hold without restriction. Condition precedent to meet the above listed requirement is the transportation in a way that an unintentionally switching on and damage by other luggage is avoided.

Therefore the **batteries must be removed*** for transportation and must be carried in strong, rigid packaging (protected from damage and from short circuit) in the cargo compartment.

The pilot-in-command must be informed about the location of the battery.

We strongly recommend obtaining the approval by corresponding Airline Company before commencement date respectively departure date.

This certificate does not apply to damaged or defect batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2022**.

Albstadt, 11th of November 2021

Alber GmbH

Hansjörg Reiner
Managing Director

Reference to battery type: **P** – Panasonic LC Series, **M** – Multipower MP Series (Jewo)

*The battery of the scala eco S34 cannot be removed. Therefore the cable to the handling unit has to be unplugged from the drive unit to avoid unintentional activation. Also the on/off switch of the drive unit has to be protected mechanically against switching on the drive unit. The approval of the corresponding Airline Company shall be obtained.