



## **Certificate in relation to air transportation by passengers with walking disabilities**

We hereby certify that the our product

### **Twion M24 / T24**

including the Lithium ion batteries fulfil the requirements of

### **IATA chapter 2.3.2.4**

and therefore may be carried as checked baggage without restriction.

The lithium-ion batteries of the twion M24 /T24 meet the requirements of the tests in the UN Manual of Tests and Criteria" Part III subsection 38.3.

The rated energy of the batteries in each wheel (83,25 Wh) is below the IATA limit of 160 Wh.

Condition precedent to meet the above listed requirement is the transportation of the twion wheels in a way that an unintentionally switching on and damage by other luggage is avoided. We recommend using the dedicated Alber transport axle.

**We strongly recommend contacting the airline sufficiently in advance to request approval for transport.**

This certificate does not apply to damaged or defective batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31<sup>st</sup> December 2021**.

Albstadt, 20<sup>nd</sup> of November 2020

Alber GmbH

A handwritten signature in black ink, appearing to be "H. Reiner", written over a light blue horizontal line.

Hansjörg Reiner  
Managing Director