



*Yes, you can.*<sup>®</sup>

# PORTFOLIO



STRENGTHEN YOUR POSITION



# CEINTURES




TOTAL LENGTH

S32 & S38	127cm (50")
M36	140cm (55")
M46	152cm (60")
L62	165cm (65")
L82 & L92	216cm (85")

ADJUSTMENT STRAP RANGE

S32 & S38	7.5cm (3")
M36	9cm (3-1/2")
M46	10cm (4")
L62	11.5cm (4-1/2")
L82 & L92	15cm (6")



WEBBING WIDTH


S32 & S38	25mm (1")
M36 & M46	38mm (1-1/2")
L62, L82, L92	50mm (2")

PADED LENGTH

S32	32cm (12-1/2")
S38	38cm (15")
M36	36cm (14")
M46	46cm (18")
L62	62cm (24-1/2")
L82	82cm (32")
L92	92cm (36")


PAD WIDTH

S32	4.5cm (1-3/4")
S38	5cm (2")
M36 & M46	5.5cm (2-1/4")
L62	7cm (2-3/4")
L82	9cm (3-1/2")
L92	10cm (4")




PUSH-BUTTON

SIZE	PART NUMBER
S32	N/A
S38	N/A
M36	HB205-M36
M46	HB205-M46
L62	HB205-L62
L82	HB205-L82
L92	HB205-L92




PLASTIC SIDE-RELEASE

SIZE	PART NUMBER
S32	HB215-S32
S38	HB215-S38
M36	HB215-M36
M46	HB215-M46
L62	HB215-L62



REHAB LATCH

SIZE	PART NUMBER
S32	N/A
S38	N/A
M36	HB225-M36
M46	HB225-M46
L62	HB225-L62



X-SMALL METAL PUSH-BUTTON

SIZE	PART NUMBER
S32	HB235-S32
S38	HB235-S38
M36	N/A
M46	N/A
L62	N/A



SWIVEL

SIZE	PART NUMBER
S32	HB275-S32
S38	HB275-S38
M36	HB275-M36
M46	HB275-M46
L62	N/A

Sub-ASIS COMPATIBLE BELT

S38	N/A
M46	LB205-M46
L62	LB205-L62

Sub-ASIS COMPATIBLE BELT

S38	N/A
M46	LB225-M46
L62	LB225-L62

Sub-ASIS COMPATIBLE BELT

S38	LB235-S38
M46	N/A
L62	N/A

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:



**-A1** + **FLAT-MOUNT™**  
END FITTING PART NUMBERS  
S (1") HKA1-R-1  
M (1-1/2") HKA1-G-1  
L (2") HKA1-B-1

**-B1** +

**-A2** + **CINCH-MOUNT™**  
END FITTING PART NUMBERS  
S (1") HKA2-R-1  
M (1-1/2") HKA2-G-1  
L (2") HKA2-B-1

**-B2**


description			
<b>HIP BELT, 2-POINT, CENTER-PULL, TECHNICAL SPECS</b>			
drawing no.	date originated	drawn by	rev.
<b>D15-0210-30</b>	2/10/2015	D. HAGAN	<b>D</b>



PLASTIC  
SIDE-RELEASE

SIZE	PART NUMBER
S32	HB216-S32
S38	HB216-S38
M36	HB216-M36
M46	HB216-M46
L62	HB216-L62

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:

<b>-A1</b>		<b>FLAT-MOUNT™</b> END FITTING PART NUMBERS S (1") HKA1-R-1 M (1-1/2") HKA1-G-1 L (2") HKA1-B-1
<b>-B1</b>		
<b>-A2</b>		<b>CINCH-MOUNT™</b> END FITTING PART NUMBERS S (1") HKA2-R-1 M (1-1/2") HKA2-G-1 L (2") HKA2-B-1
<b>-B2</b>		

note

description

**HIP BELT, 2-POINT, DUAL-PULL, TECHNICAL SPECS**

drawing no.

**D15-0210-32**

date originated


2/10/2015

drawn by


D. HAGAN

rev.

**C**




TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:



**PUSH-BUTTON**


SIZE	PART NUMBER
S38	N/A
M40	HB209-M40
M46	HB209-M46

Sub-ASIS COMPATIBLE BELT	
M46	LB209-M46



**PLASTIC SIDE-RELEASE**

SIZE	PART NUMBER
S38	HB219-S38
M40	HB219-M40
M46	HB219-M46



**REHAB LATCH**

SIZE	PART NUMBER
S38	N/A
M40	HB229-M40
M46	HB229-M46

Sub-ASIS COMPATIBLE BELT	
M46	LB229-M46



**X-SMALL METAL PUSH-BUTTON**

SIZE	PART NUMBER
S38	HB239-S38
M40	N/A
M46	N/A

**-A1** + **FLAT-MOUNT™**  
END FITTING PART NUMBERS  
S (1") HKA1-R-1  
M (1-1/2") HKA1-G-1  
L (2") HKA1-B-1

**-B1** + **FLAT-MOUNT™**  
END FITTING PART NUMBERS  
S (1") HKA1-R-1  
M (1-1/2") HKA1-G-1  
L (2") HKA1-B-1

**-A2** + **CINCH-MOUNT™**  
END FITTING PART NUMBERS  
S (1") HKA2-R-1  
M (1-1/2") HKA2-G-1  
L (2") HKA2-B-1



**-B2** + **CINCH-MOUNT™**  
END FITTING PART NUMBERS  
S (1") HKA2-R-1  
M (1-1/2") HKA2-G-1  
L (2") HKA2-B-1

note

description


**HIP BELT, 2-POINT, REAR-PULL, TECHNICAL SPECS**

drawing no. <b>D16-1116-31</b>	date originated 11/16/2016	drawn by D. HAGAN	rev. <b>A</b>
-----------------------------------	-------------------------------	----------------------	------------------



**PUSH-BUTTON**

SIZE	PART NUMBER
S38	N/A
M46	EB205-M46
L62	EB205-L62



**REHAB LATCH**

SIZE	PART NUMBER
S38	N/A
M46	EB225-M46
L62	EB225-L62



**X-SMALL METAL  
PUSH-BUTTON**


SIZE	PART NUMBER
S38	EB235-S38
M46	N/A
L62	N/A



**SWIVEL**

SIZE	PART NUMBER
S38	EB275-S38
M46	EB275-M46
L62	N/A

**INCLUDED MOUNTING HARDWARE KIT (HKA11)**



note

description

**EVOFLEX® PELVIC STABILIZER, TECHNICAL SPECS**

drawing no.

D11-1006-30

date originated

10/06/2011

drawn by

D. HAGAN

rev.

**D**

TOTAL LENGTH

S32 & S38	127cm (50")
M36	140cm (55")
M46	152cm (60")
L62	165 cm (65")

PADED LENGTH

S32	32cm (12-1/2")
S38	38cm (15")
M36	36cm (14")
M46	46cm (18")
L62	62cm (24-1/2")

CLOSED CELL FOAM PAD  
S = 7mm THICK  
M, L = 11mm THICK

TONGUE  
BUCKLE

WEBBING WIDTH

S32 & S38	25mm (1")
M36 & M46	38mm (1-1/2")
L62	50mm (2")

D-RING

ADJUSTMENT STRAP RANGE


S32 & S38	7.5cm (3")
M36	9cm (3-1/2")
M46	10cm (4")
L62	11.5cm (4-1/2")

SECONDARY STRAP LENGTH

S32 & S38	46 cm (18")
M36	51 cm (20")
M46	58.5 cm (23")
L62	63.5 cm (25")


PAD WIDTH

S32	4.5cm (1-3/4")
S38	5cm (2")
M36 & M46	5.5cm (2-1/4")
L62	7cm (2-3/4")




PUSH-BUTTON

SIZE	PART NUMBER
S32	N/A
S38	N/A
M36	HB405-M36
M46	HB405-M46
L62	HB405-L62



PLASTIC  
SIDE-RELEASE

SIZE	PART NUMBER
S32	HB415-S32
S38	HB415-S38
M36	HB415-M36
M46	HB415-M46
L62	HB415-L62



REHAB LATCH

SIZE	PART NUMBER
S32	N/A
S38	N/A
M36	HB425-M36
M46	HB425-M46
L62	HB425-L62



X-SMALL METAL  
PUSH-BUTTON

SIZE	PART NUMBER
S32	HB435-S32
S38	HB435-S38
M36	N/A
M46	N/A
L62	N/A



SWIVEL

SIZE	PART NUMBER
S32	N/A
S38	HB475-S38
M36	N/A
M46	HB475-M46
L62	N/A

Sub-ASIS COMPATIBLE BELT

S38	N/A
M46	LB405-M46
L62	LB405-L62

Sub-ASIS COMPATIBLE BELT

S38	N/A
M46	LB425-M46
L62	LB425-L62

Sub-ASIS COMPATIBLE BELT

S38	LB435-S38
M46	N/A
L62	N/A

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:

**-A1** **FLAT-MOUNT™**  
END FITTING PART NUMBERS  
S (1") HKA1-R-1  
M (1-1/2") HKA1-G-1  
L (2") HKA1-B-1

**-B1**

**-A2** **CINCH-MOUNT™**  
END FITTING PART NUMBERS  
S (1") HKA2-R-1  
M (1-1/2") HKA2-G-1  
L (2") HKA2-B-1


**-B2**

description

**HIP BELT, 4-POINT, CENTER-PULL, TECHNICAL SPECS**

drawing no. **D15-0210-31** date originated **2/10/2015** drawn by **D. HAGAN** rev. **C**

# HARNAIS DE JAMBES



SIZE	PART NUMBERS	
	FLAT MOUNT	CINCH MOUNT
SMALL	LH600S-A1	LH600S-A2
MEDIUM	LH600M-A1	LH600M-A2
LARGE	LH600L-A1	LH600L-A2

**-A1** +

FRAMESAVER CLAMP & FLAT-MOUNT END FITTING 4x  
S & M - HKA1-R-1  
L - HKA1-G-1

**-A2**

FRAMESAVER CLAMP & CINCH-MOUNT END FITTING 4x  
S & M - HKA2-R-1  
L - HKA2-G-1

SHEET 1 of 2

note

description

**LEG HARNESS TECHNICAL SPECS**

drawing no. <b>D11-0425-30</b>	date originated 4/25/2011	drawn by D. HAGAN	rev. <b>C</b>
-----------------------------------	------------------------------	----------------------	------------------




**Bodypoint**<sup>®</sup>

technical bulletin

MINIMUM LENGTH

SMALL	35-1/2" (90 cm)
MEDIUM	40-1/2" (103 cm)
LARGE	45" (114.5 cm)




PAD WIDTH

SMALL	1-1/2" (4 cm)
MEDIUM	1-3/4" (4.5 cm)
LARGE	2" (5 cm)


PAD LENGTH

SMALL	9" (23 cm)
MEDIUM	13" (33 cm)
LARGE	15" (38 cm)



MAXIMUM LENGTH

SMALL	37-1/2" (95 cm)
MEDIUM	44" (112 cm)
LARGE	50" (127 cm)



SHEET 2 of 2

note

description

**LEG HARNESS TECHNICAL SPECS**

drawing no.

D11-0425-30

date originated

4/25/2011

drawn by

D. HAGAN

rev.


**C**

# HARNAIS


# STAYFLEX



**STANDARD**



**NARROW**



**WITH ZIPPER**

	A	B	C	D	E
SMALL	9" (23 cm)	5-1/2" (14 cm)	6-1/2" (16.5 cm)	1-3/4" (4.5 cm)	4-3/4" (12 cm)
MEDIUM	12" (30.5 cm)	7" (18 cm)	8" (20.5 cm)	1-3/4" (4.5 cm)	5" (13 cm)
LARGE	15" (38 cm)	8-1/2" (21.5 cm)	9" (23 cm)	2-1/4" (5.5 cm)	5-1/4" (13.5 cm)
X-LARGE	18" (46 cm)	10" (25.5 cm)	10-1/2" (26.5 cm)	2-1/4" (5.5 cm)	6-1/2" (16.5 cm)

**WITH ZIPPER**

	A	B	C	D	E
MEDIUM	12" (30.5 cm)	7" (18 cm)	8" (20.5 cm)	1-3/4" (4.5 cm)	3-1/2" (9 cm)
LARGE	15" (38 cm)	8-1/2" (21.5 cm)	9" (23 cm)	2-1/4" (5.5 cm)	3-1/2" (9 cm)
X-LARGE	18" (46 cm)	10" (25.5 cm)	10-1/2" (26.5 cm)	2-1/4" (5.5 cm)	3-3/4" (9.5 cm)

**WITHOUT ZIPPER**

	A	B	C	D	E
X-SMALL	8-1/4" (21 cm)	5" (13 cm)	5-1/4" (13.5 cm)	1-1/2" (4 cm)	3-3/4" (9.5 cm)
SMALL	8-1/4" (21 cm)	5-3/4" (14.5 cm)	6-1/2" (16.5 cm)	1-3/4" (4.5 cm)	4" (10 cm)
MEDIUM	11-1/4" (28.5 cm)	7-1/4" (18.5 cm)	8" (20.5 cm)	1-3/4" (4.5 cm)	4-1/4" (11 cm)
LARGE	14-1/4" (36 cm)	8-3/4" (22 cm)	9" (23 cm)	2-1/4" (5.5 cm)	4-1/2" (11.5 cm)
X-LARGE	17" (43 cm)	10-1/4" (26 cm)	10-1/2" (26.5 cm)	2-1/4" (5.5 cm)	5-1/2" (14 cm)

**WITHOUT ZIPPER**

	A	B	C	D	E
MEDIUM	11-1/4" (28.5 cm)	7-1/4" (18.5 cm)	8" (20.5 cm)	1-3/4" (4.5 cm)	2-3/4" (7 cm)
LARGE	14-1/4" (36 cm)	8-3/4" (22 cm)	9" (23 cm)	2-1/4" (5.5 cm)	2-3/4" (7 cm)
X-LARGE	17" (43 cm)	10-1/4" (26 cm)	10-1/2" (26.5 cm)	2-1/4" (5.5 cm)	3" (7.5 cm)

note

description

**STAYFLEX SEWN DIMENSIONS**

drawing no. <b>D09-0729-30</b>	date originated 7/29/2009	drawn by D. HAGAN	rev. <b>B</b>
-----------------------------------	------------------------------	----------------------	------------------

SIZE	PART NUMBER
X-SMALL	SH311XS

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:

**WITH GROMMETED TOP STRAP**



-B1




FLAT-MOUNT END FITTING 2x  
HKA1-R-1

-B2



CINCH-MOUNT END FITTING 2x  
HKA2-R-1

**WITH PLAIN WEBBING TOP STRAP**



-C2-B1




FLAT-MOUNT END FITTING 2x  
HKA1-R-1

-C2-B2



CINCH-MOUNT END FITTING 2x  
HKA2-R-1



SHEET 1 of 2

note

REPLACEMENT BOTTOM STRAP * FOR ALL STAYFLEX SIZES		
SH303XL	FITS BODYPOINT SWIVEL BUCKLE ALL STAYFLEX SOLD AFTER FALL 2013	62cm (24-1/2") LONG

description

STAYFLEX™ CHEST SUPPORT, FRONT-PULL, X-SMALL, TECHNICAL SPECS

drawing no.

D11-0831-31

date originated



8/31/2011

drawn by

D. HAGAN

rev.

F



SHEET 2 of 2

note

description

STAYFLEX™ CHEST SUPPORT, FRONT-PULL, X-SMALL, TECHNICAL SPECS

drawing no.

D11-0831-31

date originated

8/31/2011

drawn by

D. HAGAN

rev.

**F**

PART NUMBERS

WITH ZIPPER SIZE	STANDARD	NARROW
SMALL	SH310S	----
MEDIUM	SH310M	SH330M
LARGE	SH310L	SH330L
X-LARGE	SH310XL	SH330XL

WITHOUT ZIPPER SIZE	STANDARD	NARROW
SMALL	SH311S	----
MEDIUM	SH311M	SH331M
LARGE	SH311L	SH331L
X-LARGE	SH311XL	SH331XL

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:

**WITH GROMMETED TOP STRAP**



-B1




FLAT-MOUNT END FITTING 2x  
HKA1-R-1

-B2



CINCH-MOUNT END FITTING 2x  
HKA2-R-1

**WITH PLAIN WEBBING TOP STRAP**



-C2-B1




FLAT-MOUNT END FITTING 2x  
HKA1-R-1

-C2-B2



CINCH-MOUNT END FITTING 2x  
HKA2-R-1



SHEET 1 of 2


note

REPLACEMENT BOTTOM STRAP * FOR ALL STAYFLEX SIZES		
SH303XL	FITS BODYPOINT SWIVEL BUCKLE ALL STAYFLEX SOLD AFTER FALL 2013	62cm (24-1/2") LONG

description

**STAYFLEX™ CHEST SUPPORT, FRONT-PULL, TECHNICAL SPECS**

drawing no. <b>D11-0825-31</b>	date originated 8/25/2011	drawn by D. HAGAN	rev. <b>F</b>
-----------------------------------	------------------------------	----------------------	------------------



MAXIMUM LENGTH

SMALL	9cm (3-1/2")
MEDIUM	9cm (3-1/2")
LARGE	11.5cm (4-1/2")
X-LARGE	11.5cm (4-1/2")

MINIMUM LENGTH  
1cm (1/2")  
ALL SIZES

BOTTOM STRAP LENGTH\*

SMALL	47cm (18-1/2")
MEDIUM	52cm (20")
LARGE	57cm (22-1/2")
X-LARGE	62cm (24-1/2")

CENTER HEIGHT

SMALL	13.5cm (5-1/4")
MEDIUM	17cm (6-3/4")
LARGE	20.5cm (8-1/4")
X-LARGE	24.5cm (9-3/4")

CENTER HEIGHT

SMALL	14.5cm (5-3/4")
MEDIUM	18.5cm (7-1/4")
LARGE	22cm (8-3/4")
X-LARGE	26cm (10-1/4")

WITH ZIPPER

WITHOUT ZIPPER

SHEET 2 of 2

note

description

**STAYFLEX™ CHEST SUPPORT, FRONT-PULL, TECHNICAL SPECS**

drawing no.

**D11-0825-31**

date originated


8/25/2011

drawn by

D. HAGAN

rev.

**F**



	CINCH-MOUNT™ END FITTING 2x HKA2-R-1
	FLAT-MOUNT™ END FITTING 2x HKA1-R-1

SIZE	PART NUMBERS	
	CINCH-MOUNT™	FLAT-MOUNT™
X-SMALL	SH321XS-B2	SH321XS-B1

SHEET 1 of 2



**note**

REPLACEMENT BOTTOM STRAP * FOR ALL STAYFLEX SIZES		
SH303XL	FITS BODYPOINT SWIVEL BUCKLE ALL STAYFLEX SOLD AFTER FALL 2013	62cm (24-1/2") LONG

**description**

**STAYFLEX™ CHEST SUPPORT, REAR-PULL, X-SMALL, TECHNICAL SPECS**

drawing no. <b>D11-0831-30</b>	date originated 8/31/2011	drawn by D. HAGAN	rev. <b>D</b>
-----------------------------------	------------------------------	----------------------	------------------



SHEET 2 of 2

note

description

STAYFLEX™ CHEST SUPPORT, REAR-PULL, X-SMALL, TECHNICAL SPECS

drawing no.

D11-0831-30

date originated


8/31/2011

drawn by

D. HAGAN

rev.

**D**



SIZE	PART NUMBERS			
	STANDARD CINCH-MOUNT™	FLAT-MOUNT™	NARROW CINCH-MOUNT™	FLAT-MOUNT™
SMALL	SH320S-B2	SH320S-B1	----	----
MEDIUM	SH320M-B2	SH320M-B1	SH340M-B2	SH340M-B1
LARGE	SH320L-B2	SH320L-B1	SH340L-B2	SH340L-B1
X-LARGE	SH320XL-B2	SH320XL-B1	SH340XL-B2	SH340XL-B1

SIZE	PART NUMBERS			
	STANDARD CINCH-MOUNT™	FLAT-MOUNT™	NARROW CINCH-MOUNT™	FLAT-MOUNT™
SMALL	SH321S-B2	SH321S-B1	----	----
MEDIUM	SH321M-B2	SH321M-B1	SH341M-B2	SH341M-B1
LARGE	SH321L-B2	SH321L-B1	SH341L-B2	SH341L-B1
X-LARGE	SH321XL-B2	SH321XL-B1	SH341XL-B2	SH341XL-B1

SHEET 1 of 2

note

REPLACEMENT BOTTOM STRAP * FOR ALL STAYFLEX SIZES		
SH303XL	FITS BODYPOINT SWIVEL BUCKLE ALL STAYFLEX SOLD AFTER FALL 2013	62cm (24-1/2") LONG

description  
**STAYFLEX™ CHEST SUPPORT, REAR-PULL, TECHNICAL SPECS**


drawing no. <b>D11-0825-30</b>	date originated 8/25/2011	drawn by <b>D. HAGAN</b>	rev. <b>E</b>
-----------------------------------	------------------------------	-----------------------------	------------------

DIRECT MOUNTING LENGTH

SMALL	40.5cm (16")
MEDIUM	40.5cm (16")
LARGE	45.5cm (18")
X-LARGE	45.5cm (18")

MAXIMUM LENGTH

SMALL	32cm (12-1/2")
MEDIUM	32cm (12-1/2")
LARGE	37cm (14-1/2")
X-LARGE	37cm (14-1/2")



ZIPPER LENGTH

SMALL	11.5cm (4-1/2")
MEDIUM	15cm (6")
LARGE	19cm (7-1/2")
X-LARGE	23cm (9")

CENTER LENGTH

SMALL	14.5cm (5-3/4")
MEDIUM	18.5cm (7-1/4")
LARGE	22cm (8-3/4")
X-LARGE	26cm (10-1/4")

BOTTOM STRAP LENGTH \*

SMALL	47cm (18-1/2")
MEDIUM	52cm (20")
LARGE	57cm (22-1/2")
X-LARGE	62cm (24-1/2")

WITH ZIPPER

WITHOUT ZIPPER

note

description

**STAYFLEX™ CHEST SUPPORT, REAR-PULL, TECHNICAL SPECS**

drawing no.

**D11-0825-30**

date originated


8/25/2011

drawn by

D. HAGAN

rev.

**E**




**PART NUMBERS**

	<b>SH351XS-B5</b> 4 - PLASTIC FLAT-MOUNT END FITTINGS HKA31-R-1
	<b>SH351XS-B6</b> 4 - PLASTIC CINCH-MOUNT END FITTINGS HKA32-R-1

**note**

<b>REPLACEMENT BOTTOM STRAP *</b> FOR ALL STAYFLEX SIZES		
SH303XL	FITS BODYPOINT SWIVEL BUCKLE ALL STAYFLEX SOLD AFTER FALL 2013	62cm (24-1/2") LONG

description <b>STAYFLEX™ CHEST SUPPORT, MD SERIES, X-SMALL, TECHNICAL SPECS</b>			
drawing no. <b>D19-0412-30</b>	date originated 4/12/2019	drawn by <b>D. DOUGAL</b>	rev. <b>A</b>



SHEET 2 of 2

note

description

STAYFLEX™ CHEST SUPPORT, MD SERIES, X-SMALL, TECHNICAL SPECS

drawing no.

D19-0412-30

date originated


4/12/2019

drawn by

D. DOUGAL

rev.

**A**



TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:

STAYFLEX WITH ZIPPER PART NUMBERS		
SIZE	STANDARD	NARROW
SMALL	SH350S	----
MEDIUM	SH350M	SH360M
LARGE	SH350L	SH360L
X-LARGE	SH350XL	SH360XL

-B5		4 - PLASTIC FLAT-MOUNT END FITTINGS HKA31-R-1
-B6		4 - PLASTIC CINCH-MOUNT END FITTINGS HKA32-R-1

STAYFLEX WITHOUT ZIPPER PART NUMBERS		
SIZE	STANDARD	NARROW
SMALL	SH351S	----
MEDIUM	SH351M	SH361M
LARGE	SH351L	SH361L
X-LARGE	SH351XL	SH361XL

SHEET 1 of 2




note

REPLACEMENT BOTTOM STRAP * FOR ALL STAYFLEX SIZES		
SH303XL	FITS BODYPOINT SWIVEL BUCKLE ALL STAYFLEX SOLD AFTER FALL 2013	62cm (24-1/2") LONG

description

**STAYFLEX™ CHEST SUPPORT, MD SERIES, TECHNICAL SPECS**

drawing no. <b>D19-0411-30</b>	date originated 4/11/2019	drawn by <b>D. DOUGAL</b>	rev. <b>A</b>
-----------------------------------	------------------------------	------------------------------	------------------




note

description  
**STAYFLEX™ CHEST SUPPORT, MD SERIES, TECHNICAL SPECS**

drawing no. <b>D19-0411-30</b>	date originated 4/11/2019	drawn by <b>D. DOUGAL</b>	rev. <b>A</b>
-----------------------------------	------------------------------	------------------------------	------------------


# HARNAIS PIVOTFIT



PART NUMBERS		
SIZE	STANDARD	DYNAMIC
XX-SMALL	SH250XXS	SH270XXS
X-SMALL	SH250XS	SH270XS
SMALL	SH250S	SH270S
MEDIUM	SH250M	SH270M
LARGE	SH250L	SH270L
X-LARGE	SH250XL	SH270XL

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:

**WITH GROMMETED TOP STRAP**



**-B1**




FLAT-MOUNT END FITTING 2x  
HKA1-R-1

**-B2**



CINCH-MOUNT END FITTING 2x  
HKA2-R-1

**WITH PLAIN WEBBING TOP STRAP**



**-C2-B1**



FLAT-MOUNT END FITTING 2x  
HKA1-R-1

**-C2-B2**



CINCH-MOUNT END FITTING 2x  
HKA2-R-1

SHEET 1 of 2

note

description

**PIVOTFIT HARNESS, FRONT-PULL, TECHNICAL SPECS**

drawing no.

**D13-0417-30**

date originated

4/17/2013

drawn by

D. HAGAN

rev.

**C**



MAXIMUM LENGTH

XX-SMALL	9cm (3-1/2")
X-SMALL	
SMALL	11.5cm (4-1/2")
MEDIUM	
LARGE	14cm (5-1/2")
X-LARGE	

MINIMUM LENGTH

1cm (1/2")  
ALL SIZES

PAD LENGTH

XX-SMALL	24.5cm (9-5/8")
X-SMALL	28cm (11")
SMALL	32.5cm (12-3/4")
MEDIUM	37cm (14-5/8")
LARGE	42.5cm (16-3/4")
X-LARGE	49cm (19-1/4")

LOWER STRAP LENGTH

XX-SMALL	40.5cm (16")
X-SMALL	
SMALL	51cm (20")
MEDIUM	
LARGE	61cm (24")
X-LARGE	


COMBINED PAD WIDTH RANGE

XX-SMALL	11.5 - 19cm (4-1/2 - 7-1/2")
X-SMALL	12.5 - 20.5cm (5 - 8")
SMALL	14 - 24cm (5-1/2 - 9-1/2")
MEDIUM	15 - 25.5cm (6 - 10")
LARGE	17 - 30cm (6-3/4 - 11-3/4")
X-LARGE	19 - 32cm (7-1/2 - 12-1/2")


AVAILABLE TOP STRAPS

ALL SIZES

GROMMETTED TOP STRAP  
SH025



PLAIN WEBBING TOP STRAP  
SH024



note

description

**PIVOTFIT HARNESS, FRONT-PULL, TECHNICAL SPECS**

drawing no.

**D13-0417-30**

date originated

4/17/2013

drawn by

D. HAGAN



rev.

**C**



# Bodypoint®

## technical bulletin



SIZE	PART NUMBERS		DYNAMIC	
	STANDARD FLAT-MOUNT	CINCH-MOUNT	FLAT-MOUNT	CINCH-MOUNT
XX-SMALL	SH260XXS-B1	SH260XXS-B2	SH280XXS-B1	SH280XXS-B2
X-SMALL	SH260XS-B1	SH260XS-B2	SH280XS-B1	SH280XS-B2
SMALL	SH260S-B1	SH260S-B2	SH280S-B1	SH280S-B2
MEDIUM	SH260M-B1	SH260M-B2	SH280M-B1	SH280M-B2
LARGE	SH260L-B1	SH260L-B2	SH280L-B1	SH280L-B2
X-LARGE	SH260XL-B1	SH260XL-B2	SH280XL-B1	SH280XL-B2

SHEET 1 of 2

note

description


### PIVOTFIT HARNESS, REAR-PULL, TECHNICAL SPECS

drawing no. <b>D13-0424-30</b>	date originated 4/24/2013	drawn by D. HAGAN	rev. <b>A</b>
-----------------------------------	------------------------------	----------------------	------------------



# Bodypoint®

## technical bulletin



SHEET 2 of 2

note

description

### PIVOTFIT HARNESS, REAR-PULL, TECHNICAL SPECS

drawing no.

**D13-0424-30**

date originated

4/24/2013

drawn by

D. HAGAN


rev.

**A**



# Bodypoint®

technical bulletin



PART NUMBERS		
SIZE	STANDARD	DYNAMIC
XX-SMALL	SH290XXS	SH295XXS
X-SMALL	SH290XS	SH295XS
SMALL	SH290S	SH295S
MEDIUM	SH290M	SH295M
LARGE	SH290L	SH295L
X-LARGE	SH290XL	SH295XL

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:

<p><b>-B5</b></p> 	<p>4 - PLASTIC FLAT-MOUNT END FITTINGS HKA31-R-1</p>
<p><b>-B6</b></p> 	<p>4 - PLASTIC CINCH-MOUNT END FITTINGS HKA32-R-1</p>

SHEET 1 of 2

note

description

## PIVOTFIT HARNESS, MD SERIES, TECHNICAL SPECS

drawing no.

**D19-0410-30**

date originated

4/10/2019

drawn by

D. DOUGAL

rev.

**A**



FRONT-PULL ORIENTATION

REAR-PULL ORIENTATION

PAD LENGTH


XX-SMALL	24.5cm (9-5/8")
X-SMALL	28cm (11")
SMALL	32.5cm (12-3/4")
MEDIUM	37cm (14-5/8")
LARGE	42.5cm (16-3/4")
X-LARGE	49cm (19-1/4")

LOWER STRAP LENGTH

XX-SMALL	44.5cm (17-1/2")
X-SMALL	
SMALL	54.5cm (21-1/2")
MEDIUM	
LARGE	65cm (25-1/2")
X-LARGE	

COMBINED PAD WIDTH RANGE

XX-SMALL	11.5 - 19cm (4-1/2 - 7-1/2")
X-SMALL	12.5 - 20.5cm (5 - 8")
SMALL	14 - 24cm (5-1/2 - 9-1/2")
MEDIUM	15 - 25.5cm (6 - 10")
LARGE	17 - 30cm (6-3/4 - 11-3/4")
X-LARGE	19 - 32cm (7-1/2 - 12-1/2")



SHEET 2 of 2

note

description

**PIVOTFIT HARNESS, MD SERIES, TECHNICAL SPECS**

drawing no.

**D19-0410-30**

date originated

4/10/2019

drawn by

D. DOUGAL


rev.

**A**


# HARNAIS TRIMELINE




SH025-A01-1 TOP STRAP, 1", FRONT-PULL



FS002 SIDE RELEASE BUCKLE, 1"



CLOSED CELL FOAM PAD  
SMALL - 7mm THICK  
M & L - 11mm THICK



SIZE	PART NUMBERS	
	CINCH-MOUNT	FLAT-MOUNT
SMALL	SH230S-B2	SH230S-B1
MEDIUM	SH230M-B2	SH230M-B1
LARGE	SH230L-B2	SH230L-B1

SHEET 1 OF 3

note

description

**TRIMLINE SHOULDER HARNESS, FRONT-PULL, TECHNICAL SPECS**

drawing no.

**D11-0824-31**

date originated




8/24/2011

drawn by

D. HAGAN

rev.

**B**



note

description

**TRIMLINE SHOULDER HARNESS, FRONT-PULL, TECHNICAL SPECS**

drawing no.

**D11-0824-31**

date originated



8/24/2011

drawn by

D. HAGAN

rev.

**B**



note

description

**TRIMLINE SHOULDER HARNESS, FRONT-PULL, TECHNICAL SPECS**

drawing no.

**D11-0824-31**

date originated

**8/24/2011**

drawn by

**D. HAGAN**


rev.

**B**



# Bodypoint®

technical bulletin



SIZE	TRIMLINE PART NUMBER
SMALL	SH235S
MEDIUM	SH235M
LARGE	SH235L

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:

-B5		4 - PLASTIC FLAT-MOUNT END FITTINGS HKA31-R-1
-B6		4 - PLASTIC CINCH-MOUNT END FITTINGS HKA32-R-1

SHEET 1 of 2

note

description

**TRIMLINE HARNESS, MD SERIES, TECHNICAL SPECS**

drawing no.

**D19-0412-31**

date originated

4/12/2019

drawn by

D. DOUGAL

rev.

**A**



FRONT-PULL ORIENTATION

REAR-PULL ORIENTATION

PAD LENGTH


SMALL	23.5cm (9-1/4")
MEDIUM	29cm (11-1/2")
LARGE	34cm (13-1/4")

PAD WIDTH

SMALL	4.5cm (1-3/4")
MEDIUM	5cm (2")
LARGE	5.5cm (2-1/4")

LOWER STRAP LENGTH

SMALL	54.5cm (21-1/2")
MEDIUM	54.5cm (21-1/2")
LARGE	65cm (25-1/2")



SHEET 2 of 2

note

description

**TRIMLINE HARNESS, MD SERIES, TECHNICAL SPECS**

drawing no.

**D19-0412-31**

date originated

4/12/2019

drawn by

D. DOUGAL

rev.


**A**



# MONOFLEX



# MONOFLEX, UNDERARM RELEASE, 1-PIECE



SIZE	W/ FLAT-MOUNT	W/ CINCH-MOUNT	W/ GROMMET STRAP
X-SMALL	SH120XS-B1	SH120XS-B2	SH120XS-B3
SMALL	SH120S-B1	SH120S-B2	SH120S-B3
MEDIUM	SH120M-B1	SH120M-B2	SH120M-B3
LARGE	SH120L-B1	SH120L-B2	SH120L-B3
X-LARGE	SH120XL-B1	SH120XL-B2	SH120XL-B3

**-B1**  **FLAT-MOUNT END FITTING 2x**  
XS & S - 1" HKA1-R-1  
M, L & XL - 1-1/2" HKA1-G-1

**-B2**  **CINCH-MOUNT END FITTING 2x**  
XS & S - 1" HKA2-R-1  
M, L & XL - 1-1/2" HKA2-G-1

note

description  
**MONOFLEX CHEST SUPPORTS, 1-PIECE & 2-PIECE**

drawing no. <b>D15-0520-30</b>	date originated 5/20/2015	drawn by D.HAGAN	rev. <b>D</b>
-----------------------------------	------------------------------	---------------------	------------------

# MONOFLEX, UNDERARM RELEASE, 1-PIECE


**MAXIMUM MOUNTING LENGTH \***

X-SMALL	68cm (26-3/4")
SMALL	76cm (30")
MEDIUM	98.5cm (38-3/4")
LARGE	106.5cm (42")
X-LARGE	114.5cm (45")


**MINIMUM MOUNTING LENGTH \***

X-SMALL	55cm (21-3/4")
SMALL	63.5cm (25")
MEDIUM	73cm (28-3/4")
LARGE	81cm (32")
X-LARGE	89cm (35")

GROMMET MOUNTING STRAP



PLAIN WEBBING MOUNTING STRAP



**MINIMUM MOUNTING LENGTH \***

X-SMALL	45cm (17-3/4")
SMALL	53.5cm (21")
MEDIUM	65.5cm (25.75")
LARGE	73.5cm (29")
X-LARGE	81cm (32")

**MAXIMUM MOUNTING LENGTH \***

X-SMALL	88cm (34-3/4")
SMALL	96.5cm (38")
MEDIUM	131.5cm (51-3/4")
LARGE	140cm (55")
X-LARGE	147cm (58")

\* NOTE: MOUNTING LENGTHS MEASURED WITH 5cm (2") ADJUSTMENT TAKEN UP AT EACH BUCKLE

note

description			
<b>MONOFLEX CHEST SUPPORTS, 1-PIECE &amp; 2-PIECE</b>			
drawing no.	date originated	drawn by	rev.
<b>D15-0520-30</b>	5/20/2015	D.HAGAN	<b>D</b>


# MONOFLEX, CENTER RELEASE, 2-PIECE

MAXIMUM MOUNTING LENGTH\*

X-SMALL	68.5cm (27")
SMALL	76cm (30")
MEDIUM	101.5cm (40")
LARGE	109cm (43")
X-LARGE	117cm (46")

MINIMUM MOUNTING LENGTH\*

X-SMALL	55cm (22")
SMALL	63.5cm (25")
MEDIUM	76cm (30")
LARGE	84cm (33")
X-LARGE	91.5cm (36")



PAD WIDTH

X-SMALL	5.5cm (2-1/4")
SMALL	6.5cm (2-1/2")
MEDIUM	8cm (3-1/4")
LARGE	9cm (3-1/2")
X-LARGE	9.5cm (3-3/4")

PADDED LENGTH

X-SMALL	31cm (12-1/4")
SMALL	39.5cm (15-1/2")
MEDIUM	46.5cm (18-1/4")
LARGE	54.5cm (21-1/2")
X-LARGE	62cm (24-1/2")

ADJUSTMENT STRAP RANGE

X-SMALL	9.5cm (3-3/4")
SMALL	
MEDIUM	11cm (4-1/4")
LARGE	
X-LARGE	

\* NOTE: MOUNTING LENGTHS MEASURED WITH 5cm (2") ADJUSTMENT TAKEN UP AT BUCKLE

SIZE	W/ FLAT-MOUNT	W/ CINCH-MOUNT
X-SMALL	SH122XS-B1	SH122XS-B2
SMALL	SH122S-B1	SH122S-B2
MEDIUM	SH122M-B1	SH122M-B2
LARGE	SH122L-B1	SH122L-B2
X-LARGE	SH122XL-B1	SH122XL-B2

**-B1**  FLAT-MOUNT  
END FITTING 2x  
XS & S - 1" HKA1-R-1  
M, L & XL - 1-1/2" HKA1-G-1

**-B2**  CINCH-MOUNT  
END FITTING 2x  
XS & S - 1" HKA2-R-1  
M, L & XL - 1-1/2" HKA2-G-1

note

description

**MONOFLEX CHEST SUPPORTS, 1-PIECE & 2-PIECE**

drawing no.

**D15-0520-30**

date originated

5/20/2015


drawn by

D.HAGAN

rev.

**D**

# CEINTURES UNIVERSELLES



**SP110S & SP111S**

note

description

**UNIVERSAL ELASTIC STRAP & EXTENSION STRAP, TECHNICAL SPECS**

drawing no.

**D13-0715-30**

date originated


7/15/2013

drawn by

D. HAGAN

rev.

**A**



**SP110L & SP111L**

note

description

**UNIVERSAL ELASTIC STRAP & EXTENSION STRAP, TECHNICAL SPECS**

drawing no.

**D13-0715-30**

date originated


7/15/2013

drawn by

D. HAGAN

rev.

**A**



note

description

**UNIVERSAL ELASTIC STRAP & EXTENSION STRAP, TECHNICAL SPECS**

drawing no.

**D13-0715-30**

date originated

7/15/2013


drawn by

D. HAGAN

rev.


**A**

# MAINTIENS DE CHEVILLES



**TABLE 1**

SIZE	PART NUMBER	SIZING GUIDE ANKLE CIRCUMFERENCE
X-SMALL	FT240XS	14 - 17cm (5-1/2 - 6-3/4")
SMALL	FT240S	16.5 - 20cm (6-1/2 - 8")
MEDIUM	FT240M	19 - 25.5cm (7-1/2 - 9-1/4")
LARGE	FT240L	23 - 28cm (9 - 11")
X-LARGE	FT240XL	26.5 - 33cm (10-1/2 - 13")



note

description

**ANKLE HUGGERS TECHNICAL SPECS**

drawing no. <b>D18-0206-30</b>	date originated 2/6/2018	drawn by D HAGAN	rev. <b>A</b>
-----------------------------------	-----------------------------	---------------------	------------------




# AEROMESH






**FRONT SIDE (AEROMESH)**



**BACK SIDE**



PART NUMBER	SIZE
SP102S	SMALL
SP102M	MEDIUM
SP102L	LARGE

SHEET 1 of 2

note

description

**AEROMESH CALF STRAP, TECHNICAL SPECS**

drawing no.

D11-0831-32

date originated


8/31/2011

drawn by

D. HAGAN

rev.

**A**



SHEET 2 of 2

note

description

**AEROMESH CALF STRAP, TECHNICAL SPECS**

drawing no.

D11-0831-32

date originated

8/31/2011

drawn by


D. HAGAN

rev.

**A**




**FRONT SIDE (AEROMESH)**



9mm AEROMESH FABRIC FACE PAD

**BACK SIDE**



VELCRO HOOK PATCH

VELCRO-COMPATIBLE LOOP FABRIC

PART NUMBER	SIZE
SP103S	SMALL
SP103M	MEDIUM
SP103L	LARGE

SHEET 1 of 2

note

description


**AEROMESH CALF PANEL, TECHNICAL SPECS**

drawing no.  
D11-0831-33

date originated  
8/31/2011

drawn by  
D. HAGAN

rev.  
**A**



SHEET 2 of 2

note

description

**AEROMESH CALF PANEL, TECHNICAL SPECS**

drawing no.

D11-0831-33

date originated

8/31/2011

drawn by

D. HAGAN

rev.

**A**



Bodypoint<sup>®</sup>  
STRENGTHEN YOUR POSITION



*Yes, you can.*<sup>®</sup>