



Product Service

CERTIFICATE

No. Z1 08 08 12941 106

Holder of Certificate: Invacare Deutschland GmbHKleiststrasse 49
32457 Porta Westfalica
GERMANY**Production Facility(ies):**

54737

Certification Mark:**Product:**

Scooters

Model(s):

Comet

Parameters:

model:	Comet 4-wheel / 3-wheel
class of use:	C
total height:	1255 mm
total length:	3-wheel: 1450 mm, 4-wheel: 1475 mm
total width:	660 mm - 725 mm
max. obstacle height:	100 mm
max. safe slope:	10°
max. speed:	10 km/h / 12.8 km/h / 15 km/h
max. load:	160 kg
backrest:	adjustable: 95°... 145°
driving motor:	model: 511100-62800 & 511100-71500
power electronics:	Dynamic Controls Ltd., Rhino DS112K
batteries:	VRLA, 2x (12 V, 63,8 Ah (5hr)), type: TEV12750
battery charger:	CTE Corporation, model: BAT-GCO812, 24 V DC, 8A (not part of testing)

Tested according to: DIN EN 12184:1999

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. The certification mark must not be altered in any way. See also notes overleaf.

Test report no.: 71340110**Date,** 2008-09-24

Page 1 of 1

