



Largeurs des dossiers Matrix selon la largeur d'assise choisie sur küschall

Largeurs d'assise	KSL	Advance	R33	K-Series	K-Series attract	Champion	Compact	Compact attract	Ultra-Light
360 mm	360 mm	360 mm	360 mm	360 mm	360 mm	360 mm	360 mm	360 mm	360 mm
380 mm	380 mm	380 mm	380 mm	380 mm	380 mm	380 mm	380 mm	380 mm	380 mm
400 mm	380 mm	380 mm	380 mm	380 mm	380 mm	380 mm	380 mm	380 mm	380 mm
420 mm	410 mm	410 mm	410 mm	410 mm	410 mm	410 mm	410 mm	410 mm	410 mm
440 mm	430 mm	430 mm	430 mm	430 mm	430 mm	430 mm	430 mm	430 mm	430 mm
460 mm		460 mm		460 mm	460 mm	460 mm	460 mm	460 mm	460 mm
480 mm				480 mm	480 mm	480 mm	480 mm	480 mm	480 mm
500 mm							480 mm	480 mm	480 mm