

Certificate in terms of transportation

We hereby confirm, that the lead-acid batteries used in our products

e-fix (E12, E15, E19, E20, E25, E26)
viamobil (V10, V15)
viamobil eco (V14)
viaplus (V12)
scalamobil (S10, S20, S24, S25, S27, S28, S30, S31, S34,
S35, S36, S38, S39)
Adventure (A10)
Quix (Q10)

are classified as

UN Nr. 2800 „Batteries, gel, nonspillable, electric storage“

according to IATA (International Air Transportation Association).

According to special provisions A67, A164 and SP872 of IATA Dangerous Goods Regulations, they are NOT RESTRICTED.

They are also considered as „nonspillable battery“ by DOT (U.S. Department of Transportation) by passing the vibration test and the pressure differential test as required in 49 CFR 173.159(d).

No electrolyte will flow from cracked case at a temperature of 55°C.

Albstadt, 2015-02-24
Alber GmbH



Ralf Ledda
Managing Director