



Product Service

CERTIFICATE

No. Z1 16 10 39801 008

Holder of Certificate: Invacare Portugal LdaRua Estrada Velha, 949
4465-784 Leça do Balio
PORTUGAL**Production Facility(ies):** 97389**Certification Mark:****Product:** General Medical Devices
Ceiling Hoist**Model(s):** Robin / Robin Mover**Parameters:**
Rated supply voltage: 24 VDC (internally powered)
Rated input power: 0.9 A
Type of applied part: B
Degree of protection: IPX4
Max. user weight: 200 kg
Mode of operation: 10 % (2 min. / 18 min.)
Lift Range: 2500 mm**Tested according to:** IEC 60601-1(ed.3);am1
IEC 60601-1-11(ed.2)
IEC 60601-1-6(ed.3);am1
IEC 62366(ed.1);am1
ISO 10535:2006

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.: 713067017**Valid until:** 2021-10-11**Date,** 2016-10-12 (Torsten Zimmer)

Page 1 of 1

